


Department of Biochemistry
Government medical college Surat
Student Journal for Practical Biochemistry



Certificate

This is to certify that _____
student of 1st MBBS, Roll No _____, Year of
Admission _____, has completed training in
practical biochemistry at Department of Biochemistry
Government Medical College, Surat.

Tutor
Department of Biochemistry
Government Medical College
Surat

Professor and Head
Department of Biochemistry
Government Medical College
Surat

INDEX

No	Exercise	Date	Page No	Sign
1	Introduction to Practical Biochemistry		04	
2	Chemistry of Carbohydrates		06	
3	Chemistry of Proteins and Amino acids		17	
4	Chemistry of lipids		27	
5	Physiological urine		30	
6	Pathological urine		37	
7	Estimation of acid output by stomach		44	
8	Secretion and buffering of acids by kidney		51	
9	Colorimetry		56	
10	Estimation of Serum Creatinine		61	
11	Estimation of Plasma Glucose		64	
12	Estimation of Serum Cholesterol		66	
13	Estimation of Serum Bilirubin		68	
14	Estimation of Serum Total Protein		71	
15	Estimation of Serum Albumin		73	
16	Estimation of Cerebro Spinal fluid protein		76	
17	Estimation of serum uric acid		78	
18	Electrophoresis		80	
19	Chromatography		83	
20	Clinical Cases		85	
21	Medical Biochemistry – Syllabus		94	
22	Subject distribution and Paper Style		96	

1. Introduction to Practical Biochemistry

Practical Biochemistry serves several purposes to medical students.

Practical Biochemistry augments concepts learnt in classroom. e.g.

Study of properties of basic biomolecules. e.g Carbohydrate, Proteins

Study of biochemical investigative tools e.g colorimetry, chromatography, electrophoresis

Study of patient case histories in light of its laboratory investigations is fundamental to understanding medical aspects of biochemistry.

It prepares the student for possible use of the practical techniques in clinical practice. e.g.

Many bedside biochemistry diagnostic technologies are used by physicians themselves. Such technologies are called point-of-care-technologies (**POCT**). They are based on many simple concepts studied in practical biochemistry.

Many biochemistry diagnostic technologies are used by patients themselves. Such home monitoring by patients require support from their physicians. Practical biochemistry help medical students for supporting their patient's for such support.

Hazards in Clinical Biochemistry laboratory

Hazards arises from three main basic sources

1. From dangerous chemicals
2. From infected specimen sent for analysis
3. From faulty apparatus & instruments

These are further increased by carelessness, untidiness, faulty hygiene, conduct of staff, unsatisfactory working condition.

Wide variety of articles are used like conical flasks, volumetric flasks, tube, measuring cylinders, pipettes, reagent bottles.

Disposal of laboratory Waste

There are guidelines to dispose waste. It is recommended that waste should be segregated at the point of generation & disposed in bags with correct colour coding.

YELLOW BAGS	RED BAGS	BLUE BAGS	BLACK CARBOY
Infectious waste, bandages, gauze, cotton or any other objects in contact with body fluids, human body parts, placenta etc.	Plastic waste such as catheters, injection syringes, tubings, iv bottles	All types of glass bottles and broken glass articles, outdated & discarded medicines	Needles without syringes, blades, sharps and all metal articles.

Questions:

Describe any laboratory accident you or your schoolmate has suffered in your school days. How will/was it be first-aid? How will you prevent it?

Give list of Biochemistry POCT and home-monitoring technologies. Explain each of them.

Mention five more points to be noted for laboratory safety.

2. Chemistry of Carbohydrates

Test solution

Glucose solution(400mg/dl) : Dissolve 4 gm of glucose powder in 1000 ml water

Starch solution(1%): Add 10 gm of starch powder in 100 ml of water .Boil till solution become clear. Make up to 1 litre.

Sucrose solution(400mg/dl) : Dissolve 4 gm of Sucrose powder in 1000 ml water

Fructose solution(400mg/dl) : Dissolve 4 gm of Fructose in 1000 ml water

Maltose solution(400mg/dl) : Dissolve 4 gm of Maltose powder in 1000 ml water

Molisch's test:


Reagent

1 % α -Naphthol: Dissolve 1 gm α -Naphthol powder in 100 ml methanol
Conc. H_2SO_4

Principle

All carbohydrates when treated with conc. sulphuric acid undergo dehydration to give fufural compounds. These compounds condense with Alpha-naphthol to form colored compounds.

Molish test is given by sugars with **at least five carbons** because it involves fufurl derivatives, which are five carbon compounds.



Benedict's Test:

All Reducing sugars give positive benedict's test.Reducing sugars have a free aldehyde or keto group.

Reagent

Benedict's Reagent:One liter of Benedict's solution contains ,
 173 grams -----> sodium citrate,
 100 grams -----> sodium carbonate
 17.3 grams ----->cupric sulphate pentahydrate.


With the help of heat,dissolve 173 gm of sodium citrate & 100 gm of sodium carbonate in 800 ml of water.Dissolve 17.3 gm cupric sulphate pentahydrate in 100 ml of water in different container.

Pour cupric sulfate solution in carbonate- citrate solution with constant stirring& make upto 1000ml.

Role of ingradient of benedict's solution:

- 1.Sodium citrate:Holding of cupric oxide in alkaline solution
- 2.Sodium carbonate:provide alkaline pH
- 3.cupric sulphate pentahydrate:Reducing Agent

Principle




The principle of Benedict's test is that when reducing sugars are heated in the presence of an alkali(pH 10.6), they get converted to powerful reducing

compounds known as enediols. Enediols reduce the cupric ions (Cu^{2+})

present in the Benedict's reagent to cuprous ions (Cu^+) which get precipitated as insoluble red copper(I) oxide.

The color of the obtained precipitate gives an idea about the quantity of sugar present in the solution, hence the test is semi-quantitative.



Carbohydrates giving positive Benedict's test:

Glucose, Fructose, Galactose,
Ribose, Glucuronic acid,
Lactose, Maltose

Note: Sucrose with no free reducing group give negative test.

Non-Carbohydrates giving positive Benedict's test:

High concentration of Uric acid , Creatinine and Ketones
Homogentisic acid (solution turns black due to black colored oxidized
homogentisic acid)
Vitamin C (even without Boiling)
Certain drugs like aspirin, cephalosporins

Starches

Starches do not react or react very poorly with Benedict's reagent, due to the relatively small number of reducing sugar moieties, which occur only at the ends of carbohydrate chains.

Different concentration of glucose gives different color of solution with Benedict's test, depending on amount of precipitate and residual cupric sulphate.

Grade	Color of Reaction Mixture	Approximate Glucose concentration
+	Green	0.5-1 gm%
++	Yellow	1-1.5 gm%
+++	Orange	1.5-2 gm%
++++	Red	>2 gm%

Benedict's test is frequently used to detect glucose in urine. Although glucose is most frequent reducing substance present in urine, in some patient positive Benedict's test may be due to non-glucose reducing substances listed above. This phenomenon may be called **false positive** result.




Following test based on glucose oxidase is positive only with glucose in urine.

Glucose oxidase test:

Reagent:

Glucose strip or liquid reagent based on GOD-POD method

Principle





Some of the dyes used are O-tolidine, tetramethylbenzidine, and potassium iodide, 4-aminophenazome + phenol .

Reagents for this test are present on a strip of paper in solid form. When the paper is wet with urine, the reagents dissolve in urine on paper and react with glucose in urine. The darkness of color can be correlated with amount of glucose present in urine.

Because **Glucose oxidase enzyme can act only on beta-D-Glucose**, other reducing substances do not give this test positive. (Exception: Galactose can react with glucose oxidase, but very slowly)

Following reaction occur when urine contain compounds reacting with H₂O₂.



Thus, compounds like Vitamin C, Aspirin utilize H₂O₂ produced in the reaction. Due to lack of H₂O₂, peroxidase can not oxidize dye. Thus, glucose may not be detected even if present, if urine contain Vitamin C or Aspirin in large amount. This phenomenon is called **false negative** result.

In neonate, **positive Benedict's test in urine, in presence of negative Glucose oxidase test**, indicate possible presence of Fructose or Galactose in urine. (But note the exception mentioned above). Fructose and galactose are found in some inborn deficiency of enzymes of their metabolic pathways.

Performed Benedict's test and Glucose oxidase strip test with following compounds and fill up the table given.

Compound	Benedict test	Glucose Oxidase	Inference
----------	---------------	-----------------	-----------

		Strip Test	
Fructose			
Vitamin C			
Glucose with Vitamin C			
Cephalosporin Drugs			


Barfoed's Test:

This test is based on the same principle as Benedict's test. But, the test medium is acidic. In acidic medium (pH 4.6) **monosaccharides** react faster than disaccharide. Monosaccharides react fast within 1-2 minutes but disaccharides take longer i.e. 7-12 minutes.

Reagent:

Barfoed's reagent: Dissolve 70 gm of cupric acetate monohydrate in 800 ml of water. Add 9 ml glacial acetic acid & make to 1000 ml with water.

Principle



Seliwanoff's Test

Seliwanoff's test is a chemical test which distinguishes between aldose and ketose sugars. This test is based on the fact that, when heated, ketoses are more rapidly dehydrated than aldoses.


Reagent

Seliwanoff's reagent: Add 50mg of Resorcinol in 66 ml of water. Add 33 ml concentrated HCL. Wear goggles. The reagent is colorless if red color develop, discard it.

Principle

Ketohexoses like fructose on treatment with HCl form 5-hydroxymethylfurfural, which on condensation with resorcinol gives a cherry red complex.

Sucrose is hydrolyzed into glucose and fructose when boiled in acidic medium of Seliwanoff's reagent. Fructose, present in hydrolysate gives positive Seliwanoff's test.



Inversion Test:

Reagent

10 % HCL


40% NaOH : dissolve 40 gm of NaOH pellet in 100ml Water

Benedict's reagent

Seliwanoff's reagent

Principle

When sucrose is boiled with conc. HCl, It is hydrolyzed into its constituent monosaccharides i.e. fructose and glucose. The hydrolyzed glucose and fructose give Benedict's test. Fructose gives seliwanoff's test.



Sucrose is dextrorotatory. The optical rotation changes from dextrorotatory to leavorotatory on hydrolysis, since fructose causes a much greater leavorotation than the dextrorotation caused by glucose. This is known as inversion. The resultant hydrolysate is called invert sugar, which is sweeter than sucrose because fructose is sweeter than sucrose.

Iodine test for starch

Reagent:


Iodine solution : Dissolve 1.27 gm Iodine and 3 gm potassium iodide crystals in 100 ml water. Dilute 1:10 in water before use.

Iodine by itself is very poorly soluble in water. One way to dissolve iodine in water is to add potassium or sodium iodine. Those salts dissolve into potassium or sodium ions and iodine ions. The iodine ion (I⁻) reacts with the free iodine (I₂) to form a triiodide ion (I₃⁻) which is soluble in water and can react with glucose chains.

Principle

Iodine binds starch to give blue colored complex.

When glucose chains are sufficiently long they coil up like springs. This coil is supported by weak links between the glucose molecules. These links break down at high temperatures and the glucose chains uncoil. When the chains are longer than about 9 glucose molecules a triiodide ion (I_3^-) fits inside the coil (Figure).The longer the glucose chains are the more iodine molecules fit into the coils and the more intense the color reaction will be.



The resulting color depends on the length of the glucose chains. Shorter chains (starting at about 9 glucose molecules in unbranched chains and up to 60 glucose molecules in branches chains) give a red color .

Amylose, which consists of very long glucose chains between occasional branch points and very large **dextrines** give a dark blue color .

while **amylopectin**, which has much more branch points and shorter glucose chains between these branch points, gives a more reddish color in the presence of iodine.

Hydrolysis Test for starch

When starch/dextrin is boiled with HCl, It is hydrolyzed into its constituent monosaccharides i.e. glucose. Glucose, thus formed, gives Benedict's test.

TEST	METHOD	OBSERVATION	INFERENCE
Molisch's Test	1ml OS + 2 drops of α -naphthol solution mix. add 2 ml. of conc. Sulphuric acid carefully through the side of the test tube without shaking.	Purple ring is formed at the junction of acid and solution.	Carbohydrate present.
Benedict's Test	5ml of Benedict's reagent + 8 drops of OS, mix Boil and cool.	Green / Yellow / Orange / Red / Brick Red precipitates seen	Reducing Group present.
Barfoed's Test	1 ml OS + 1 ml Barfoed's reagent Boil for 30 sec , Cool Excess boiling or may give false positive results.	Red colored precipitates. At the bottom of the tube.	Disaccharides absent. Monosaccharide present
Seliwanoff's Test	1 ml O.S. + 1 ml Seliwanoff's reagent. Boil and cool for 5 min	Red colored formed.	Keto sugars present e.g. Fructose
Iodine Test	1 ml OS + 2 drops of iodine solution, Mix	Blue color develops. Violet colour develops.	Starch present. Dextrine present.
Inversion Test	1 ml OS + 1ml of 10% HCl. Boil for 2 mins. Cool. Make it alkaline with 5 drops of 40% NaOH. From this solution perform Benedict's test and Saliwanoff's test.	Benedict's and Saliwanoff's test are positive	Sucrose is present if OS give negative Benedict's test.
Hydrolysis test for starch/dextrin	Step-1: Perform Benedict's Test with OS. Step-2: 5 ml OS + 2 drops of conc. HCl . Boil for 2 mins. Cool. Make alkaline with 5 drops of 40% NaOH. From this solution perform	Benedict's test is negative/ weakly positive Benedict's test is positive	Starch present (weak Benedict's test with OS is due to free reducing groups at end of starch molecules.)

	Benedict's test		
Glucose oxidase test (on strip or with liquid reagents)	Method for the test will be provided in the laboratory	Observation will be explained in the laboratory	Glucose present in the solution

What you will do:

Perform tests mentioned in above table with various carbohydrates given to you. Note down your observation and inference in tables as shown below.

TEST	OBSERVATION	INFERENCE
Molisch's Test		
Benedict's Test		
Barfoed's Test		
Seliwanoff's Test		

Iodine Test		
Inversion Test		
Hydrolysis test for starch/dextrin		
Glucose oxidase test (on strip or with liquid reagents)		

Questions:

Explain biochemical reason why Sucrose gives negative Benedict's test.

Why the hydrolysis of sucrose is called 'Inversion test'?

Does alpha-D-Glucose in the solution react with Glucose Oxidase?
Explain.

3. Chemistry of Proteins and Amino acid

Proteins are made up of amino acids. Amino acids differ from each other in their side chain (-R group). The differing -R groups in different amino acids are responsible for many reactions mentioned below.

Preparation of Protein solutions:

Egg albumin solution (1:21): Mix 50 ml of egg(both white and yellow) in 1 liter of tap water. Use only for 24 hours

Gelatine solution(0.5%): Dissolve 5 gm of Gelatin powder in 50 ml of water by slight Heating & make upto 1 liter

Peptone solution(0.5%) :Dissolve 5 gm of Peptone powder in 50 ml of water by slight Heating & make upto 1 liter

Casein solution(0.5%) : Dissolve 5 gm of Casein powder in 20 ml of 40% NaOH & make upto 1 Liter with water

Biuret Test


This test is given by all peptides having at least two peptide bonds. So, it is given by all proteins.

Reagents:

10% NaOH: Take 10gm NaOH pellets and make it up to 100ml with DI water.

1% CuSO₄: 1 gm of CuSO₄ in 100 ml DI water.

Cu²⁺ - peptide complex



Principle:

Cupric ions of copper sulphate solutions in alkaline medium form coordinate complex with at least two nitrogens of the peptide bonds to form purple colored complex. Thus color intensity is proportionate to the presence of number of peptide linkages.

Minimum of 2 peptide bonds (3 amino acids) are required for binding of Cu^{2+} with peptide. single amino acids and dipeptides do not give positive test.

The name of reaction is derived from organic compound **biuret** which is formed by condensation of 2 urea molecules at high temperature.




Figure of Biuret

Biuret is formed when solid urea powder is heated in a tube. The resultant Biuret is solid at room temperature and soluble in water.

The test produces color proportionate to number of peptide bonds which can be correlated with amount of protein. Similar reagent is used for estimation of serum proteins quantitatively.

Ninhydrin Test

This test is given by all compounds having free α -Amino groups. ex: peptides, proteins, free α - Amino acid. Different Proline and hydroxyproline give yellow color in this test.

Prepare reagent:

1 % Ninhydrine solution : 1 gm of Ninhydrine powder dissolved in 100 ml DI water.

Principle:

Ninhydrine + α - Amino acid \rightarrow hydrindantin + aldehyde + CO₂ + NH₃
Hydrindantin + NH₃ + Ninhydrine \rightarrow blue colored complex

Ninhydrin oxidises an α -amino acid to an aldehyde liberating NH₃ and CO₂ and is itself reduced to hydrindantin. Hydrindantin then react with NH₃ and another molecule of ninhydrine to form a purple colored complex.

All amino acids that have a free amino group will give positive result (purple color) .While not free amino group-proline and **hydroxy-proline** (amino acids) will give a (yellow color).

Note: Many substances other than amino acids, such as amines will yield a blue color with ninhydrin, particularly if reaction is carried out on filter paper.

Xanthoproteic Test:

This test is answered by aromatic amino acids. (Tyrosine, Tryptophane)


Reagent:

Concentrated HNO₃

40 % NAOH : 40 gm NAOH in 100 ml DI water.

Principle

Concentrated nitric acid causes nitration of activated benzene ring of tyrosine and tryptophan. (Benzene ring is considered activated when additional groups are attached to it) The nitrated activated benzene is yellow in color. It turns to orange in alkaline medium.. Phenylalanine also contains benzene ring, but ring is not activated, so it does not undergo nitration. The reaction can be hastened by heating. The heat may be produced by dilution of concentrated HNO₃ with OS or may require heating.



Aldehyde Test

Reagents

1:500 Formaldehyde Reagent:

Take 1 ml of Formaldehyde solution (37-41 % W/V) and make upto 500 ml with DI water. Use only for 1 week. Old Formaldehyde may not give test.

1 % Sodium Nitrite solution :


Take 1 gm sodium nitrite powder and make upto 100 ml with DI water. Use only for 1 week. Old Sodium nitrite may not give test.

Sulphuric acid AR :

Use sulphuric acid Bottle directly for use as reagent. Use for 1 week. Old Sulphuric acid may not give test

Principle

Indole ring is present in tryptophan. Formaldehyde react with indole ring to give violet colored complexes in presence of H₂SO₄. Addition of Sodium nitrite intensify and stabilize colour.



Millon's reagent


Reagent:

Millon's reagent:

Dissolve 10 gm of mercuric sulphate(HgSO_4) +100ml DI water + 7 ml Conc. H_2SO_4

1% sodium nitrite: 1 gm in 100 ml DI water

Principle



Tyrosine has hydroxyphenyl(Phenol) group. The hydrophobic group is in the core of protein. The protein is denatured by mercuric sulphate in boiling water exposing hydroxyphenyl group. Sodium nitrite reacts with sulfuric acid to form nitrous acid. The exposed hydroxyphenyl groups react with nitrous acid. The compound formed chelates Hg^{2+} & give red colour precipitates.

Sakaguchi's Test

This test is for Guanido group Which is the R-group of arginine.

Reagent:

1%w/v α -Naphthol: Dissolve 1 gm α -Naphthol in 100 ml of methanol

10%w/v NaOH: Dissolve 10gm of NaOH & make it upto 100ml with DI water.

Alkaline hypochloride : Make 100 ml 10 % NaOH & add 8 ml 5-6 % Analytical grade Sodium hypochloride.

Principle

In an alkaline medium, alpha-Naphthol combines with guanidino group of arginine to form a complex, which is oxidized by bromine/chlorine.


Sulphur Test (Lead acetate test):

Reagent:

2% Lead acetate in 10% NaOH: add 20 gm lead acetate, 100 gm NaOH in 1 liter of water. There is no need to make exactly up to 1 liter. Above solution will be more than 1 liter in volume.

Principle:

When protein containing cysteine & cystine is boiled with strong alkali, organic sulphur(R-SH) is converted to sulphide (Na_2S). Addition of lead acetate to this solution causes precipitation of insoluble lead sulphide (PbS), which is black-gray in colour. Methionine does not give this test due to the presence of thioether linkage ($\text{H}_3\text{C-S-CH}_2\text{-R}$) which does not allow the release of sulphur in this reaction.



Heat coagulation test:

Reagent:

1% acetic acid: 1 ml acetic acid up to 100 ml with DI water.

Principle

Proteins have net zero charge at their iso-electric pH (pI). So, at pI, protein molecules have minimum repelling force. Thus proteins are easily precipitated at pI. When proteins are heated, weak bonds like hydrogen-bonds, salt bonds and van-der-wal forces are broken. Proteins are said to be denatured.

Core hydrophobic regions of denatured Albumin can form intermolecular associations and cause precipitation. Thus, in order to precipitate proteins like albumin, two conditions are required. 1) Bring albumin to its pI(5.4) by adding few drops of 1% acetic acid. 2) Heat the solution



Half & Full Saturation Test:


Reagent:

Saturated ammonium sulphate [(NH₄)₂SO₄]: Add ammonium sulphate in 500 ml DI water till it stops dissolving.

Ammonium sulphate [(NH₄)₂SO₄] power

Principle

When ammonium sulphate is added to protein solution, water concentration decreases. This removes shell of water from outer surface of protein molecules, favoring formation of hydrogen bonds among protein molecules and causing their precipitation. While proteins like globulin, gelatin and casein are precipitated in half-saturated ammonium sulphate solutions, albumin is precipitated in full-saturated ammonium sulphate solution.



Protein molecules contain both hydrophilic & hydrophobic aminoacids.

In aqueous medium, hydrophobic amino acids form protected areas while hydrophilic amino acids form hydrogen bonds with surrounding water molecules (solvation layer). When proteins are present in salt solutions (e.g. ammonium sulfate), some of the water molecules in the solvation layer are attracted by salt ions. When salt concentration gradually increases, the number of water molecules in the solvation layer gradually decreases until protein molecules coagulate forming a precipitate; this is known as “salting out”.

For example, albumin requires higher salt concentration for precipitation than casein or gelatin. Albumin particles are smaller in size & so have larger surface area, so they hold more water molecules around them. So a higher concentration of Ammonium sulphate is required. The salt concentration used is described as 'half saturation' (for

casein, gelatin, globulin) or 'full saturation' (for albumin).

PROCEDURES

TEST	METHOD	OBSERVATION	INFERENCE
BIURET TEST	<ul style="list-style-type: none"> ➤ 10% NaOH (2 ml) + 1% CuSO₄ (2 ml) ➤ divide above mixture in two parts of 2 ml ➤ part 1: add 2 ml OS ➤ part 2: add 2 ml H₂O 	Pink or Violet Colour develops in part 1. No such color develop in part 2	Two or more peptide linkages present. Protein present
XANTHO-PROTEIC TEST	<ul style="list-style-type: none"> ➤ OS (0.5 ml) + HNO₃_{con} (1 ml) Mix it. (Solution turns yellow) + 40%NaOH (1 ml) in above mixture. ➤ Solution turns orange <p>Note: Use Fresh(tightly packed) conc.HNO otherwise test come negative.</p>	Yellow-Orange colour develops.	Aromatic Amino Acids Tyrosine and Tryptophan present in protein.
NINHYDRIN TEST	<ul style="list-style-type: none"> ➤ OS (1 ml) + 1% Ninhydrine (2 drops) ➤ Mix, Boil (1 min). ➤ Cool. 	Blue or Purple colour develops.	Alpha Amino groups of proteins at N-terminal are responsible for positive test with proteins.

Aldehyde Test	<ul style="list-style-type: none"> ➤ 1 ml Protein Solution + 1 drop of 1:500 formalin. Mix. ➤ Slant the test tube and slowly add 1 ml of conc. H₂SO₄ . Mix. ➤ Add 1 drop of 1% sodium nitrite solution in Test tube. Mix. ➤ Use Fresh(tightly packed) conc.H₂SO₄ & 1:500formaline otherwise test come negative. 	Violet color is formed.	Indole group present in protein. Tryptophan present in the protein.
MILLION'S TEST	<ul style="list-style-type: none"> ➤ 0.5 ml protein sol. +50 ul sodium nitrate sol.ⁿ+100 ul Millon's reagent. mix well & Heat 	Red coloured precipitate Observed.	Hydroxyphenyl group present in protein. Tyrosine present in protein.
SAKAGUCHI TEST	1 ml Protein sol. ⁿ + 2 drops of alpha Naphthol + 1 ml Alkaline sodium Hypoochloride	Carmin Red colour observed.	Guanidino group present in protein. Arginine present in protein.
SULPHUR TEST (Lead acetate test)	<ul style="list-style-type: none"> ➤ 0.5 ml OS + 0.5 ml Lead acetate reagent ➤ Boil for 1 minute 	Black- Grey colour seen.	Sulfhydryl group (-SH) present in protein. Cysteine & Cystine present in protein
HEAT COAGULATION TEST	<p>5 ml Protein solution + 1- 2 drops of chlorophenol red.</p> <ul style="list-style-type: none"> ➤ If faint pink colour appears, add 1% acetic acid drop wise till you get faint pink colour with yellowish tinge (pH 5.4 is obtained) ➤ If yellow colour appears, add 2% Na₂CO₃ solution drop wise till you get faint pink colour with yellowish tinge (pH 5.4 is obtained) 	White precipitates seen in upper part of solution, as compared to clear lower part of solution	Albumin is precipitated when denatured at its pI~5.4

HALF SATURATION TEST	2 ml of the protein sol. ⁿ + 2 ml of saturated sol. ⁿ of (NH ₄) ₂ SO ₄ (Thus, saturated (NH ₄) ₂ SO ₄ is half diluted)	White precipitate formed.	Casein, Gelatin and Globulin are precipitated at half saturation with (NH ₄) ₂ SO ₄
FULL SATURATION TEST	5 ml. Of protein sol. ⁿ + a pinch of Ammonium Sulphate powder, Shake Repeat above steps till some undissolved (NH ₄) ₂ SO ₄ remains at the bottom of the test tube.	White precipitate formed	Albumin precipitates at full saturation with (NH ₄) ₂ SO ₄
MOLISCH'S TEST	<ul style="list-style-type: none"> ➤ 1ml OS + 2 drops of α-naphthol solution, mix ➤ Add 2 ml. of conc. Sulphuric acid carefully through the side of the test tube without shaking. 	Purple ring is formed at the junction of acid and solution.	Proteins contain Carbohydrates

What you will do:

- Perform tests mentioned in above table with various Protein Solutions given to you. Note down your observation and inference in tables as shown below.

TEST	OBSERVATION	INFERENCE
Biuret Test		
XANTHO-PROTEIC TEST		
NINHYDRIN TEST		
Aldehyde Test		

MILLION'S TEST		
SAKAGUCHI'S TEST		
SULPHUR TEST		
HEAT COAGULATION TEST		
HALF SATURATION TEST		
FULL SATURATION TEST		
MOLISCH'S TEST		

➤ Fill up the table given below.

Use: 'P' for positive test 'N' for negative test 'W' for weakly positive test

Test	Amino acids responsible for the test	Albumin	Casein	Gelatin	Peptone
Xanthoproteic test					
Ninhydrin test					
Hopkin's and Cole test					
Million's test					
Sakaguchi's test					
Lead acetate test					

➤ Mention food sources of Albumin, Casein and Gelatin.

➤ Which of the Albumin, Casein and Gelatin is nutritionally best? Explain.

➤ If by mistake Ninhydrin touches your skin while doing the ninhydrin test, skin gets bluish stain. Explain.

4. Chemistry of lipids

Lipids are heterogeneous group of compounds soluble in non-polar solvents like chloroform but not soluble in polar solvents like water.

While body is water medium, lipids of body require specialized methods for digestion, absorption and transport.

Bile salts cause emulsification of oil due to their amphipathic nature and ability to reduce surface tension. Thus making bile salts essential for digestion and absorption of lipids of food.

Lipids of blood are transported as lipoproteins. Without lipoproteins, lipids would be insoluble in plasma (93% water).

Reagent

Any oil : Ground nut oil, coconut oil

Non polar Solvent : Acetone/ Methanol

Bile salt solution : Dissolve 0.6 gm sodium deoxycholate in 100 ml DI water. Do not take tap water for making bile salt solution, precipitation occurs due to interference by calcium.


TEST	METHOD	OBSERVATION
Solubility of oil in water	➤ 0.1 ml of oil + 1 ml water, mix, for 15 sec.	Big oil drops are observed
Solubility of oil in non-polar solvent	➤ 0.1 ml of oil + 1 ml Acetone/Methanol, mix, for 15 sec.	oil droplets are not observed
Emulsification of oil in Bile salts.	<ul style="list-style-type: none"> ➤ Take 2 test tubes T1 and T2 ➤ Take 1 ml H₂O in T1 test tube. ➤ Take 1 ml Bile salt solution (Sodium deoxycholate solution) in T2 test tube. ➤ Add 0.1 ml of oil in T1 and T2. ➤ Mix T1 & T2, all together for 15 sec. against palm of your hand. 	<p>Compare size of oil drops and turbidity immediately .</p> <p>T1: Big oil drops, Clear water (Compared to T2)</p> <p>T2: Small oil drops, Turbid solution (Compared to T1)</p>

What will you do:


➤ Perform the test shown above with the oil provided. Draw table showing the tests and your observations.


TEST	OBSERVATION	INTERFERENCE
Solubility of oil in water		
Solubility of oil in non-polar solvent		
Emulsification of oil in Bile salts.		

Draw structure of bile salt




Draw structure of an oil droplet in a bile salt solution.






Draw structure of a micelle. Write its function in body.



Draw structure of a lipoprotein particle. Write its function in body.



5. Physiological Urine

Artificial Urine sample:

Ammonium sulfate	2 gm
Sodium phosphate dibasic(monobasic	2 gm
Pottasium dihydrogen phosphate	2 gm
Urea powder	2 gm
creatinine powder	2 gm
Uric acid powder	1 gm
Calcium carbonate/Calcium chloride :	1 gm.
NaCl	4 gm,
And make upto	2 liters

Physical characteristics of urine

Volume:

Normal adult excretes 800-2000 ml of urine daily

Factors affecting urine volume:

- According to quantity of fluid ingested
- Environment temperature
- Physical activity
- Loss of water in feces, via skin, in vomitus etc.

Collection of urine to measure volume:

Discard the first morning urine. Then collect urine during each micturition in a vessel up to, including the next morning urine.

Some conditions with increased urine volume:

- o Diabetes mellitus
- o Diabetes insipidus (low specific gravity of urine)
- o Diuretics drug therapy

Some conditions with decreased urine volume:

- o Dehydration
- o Renal failure

Appearance:

Normal urine is clear and transparent when freshly voided. On standing bacterial urease converts urea into CO₂ and Ammonia. Ammonia makes urine alkaline. Phosphates precipitate in alkaline urine making it turbid.

Color:

Fresh urine is amber yellow. This colour is due to urobilin.

Odour:

Fresh urine has an aromatic odor due to presence of volatile organic acids produced by body and intestinal bacteria.

Reaction:

Fresh urine is normally acidic (pH<7.0). Post-prandial urine is alkaline due to secretion of HCl in stomach, the condition known as “Alkaline Tide”.

Specific gravity:

Normal range-1.003 to 1.035 gm/ml of urine. The greater the amount of solutes per unit volume of urine, the greater the specific gravity. It is high in diabetes mellitus, while low in diabetes insipidus.

Determination of specific gravity: Wipe the urinometer by a filter paper and allow it to float in the urine contained in the cylinder. See carefully that the apparatus do not touch the sides or bottom of the cylinder, when it is at rest take the reading from lower meniscus (true surface) of urine. Note the temperature of urine. If it differs from the standard temperature written on the urinometer, add one unit (0.001) for every 3 degree rise from the standard temperature.

Chemical analysis of urine

A. Inorganic Chemical constituents


Ammonia:**Reagent:**

1% phenolphthalein : Dissolve 0.5 gm of phenolphthalein in 50 ml of methanol. Phenolphthalein is insoluble in water

2% sodium carbonate : Dissolve 10 gm of sodium carbonate in 500 ml of water

Principle:


Urinary ammonia is derived from glutamine in kidney. It is secreted as a buffer against H⁺ secreted by tubules.



On heating NH₃ evaporate, dissolve in water around a glass rod and make it alkaline. At alkaline pH phenolphthalein ions are formed which is pink coloured.

Phenolphthalein is a weak acid, which can lose H⁺ ions in solution. The phenolphthalein molecule(HIn) is colorless, and the phenolphthalein

ion(In⁻) is pink. When a base is added to the phenolphthalein, the molecule \rightleftharpoons ions equilibrium shifts to the right, leading to more ionization as H⁺ ions are removed.



For phenolphthalein: pH 8.2 = colorless; pH 10 = red

Procedure:

- Take 5ml urine in a test tube and add a drop of phenolphthalein. Add drop wise 2% sodium carbonate solution till the solution turns faint pink. Boil and hold a glass rod dipped in phenolphthalein at the mouth of the test tube. Phenolphthalein turns pink due to gaseous ammonia.

Chloride:

Reagent:

Concentrated HNO₃

3% AgNO₃ : Dissolve 15 gm of AgNO₃ in 500 ml of water.

Principle:



When acidified urine reacts with silver nitrate, a white precipitate of silver chloride is formed.

Procedure:

- [3 ml of urine] + [1.0ml concentrated HNO₃] + [1.0 ml 3% AgNO₃]
Curdy white precipitate of AgCl is formed.

(Concentrated HNO₃ is added to prevent precipitation of urate and acid phosphates by AgNO₃)

Calcium:

Reagent:

Saturated ammonium oxalate solution: Dissolve ammonium oxalate powder in 500 ml of water till it become undissolved.

Principle:

Calcium precipitated as insoluble calcium oxalate with ammonium oxalate



Procedure:

➤ Sulkowitch Test: To 5 ml urine and add 3 ml saturated ammonium oxalate solution.

Calcium precipitated as insoluble calcium oxalate is observed as turbidity.


Phosphorus:

Reagent:

Concentrated HNO₃

5% Ammonium Molybdate : Dissolve 5 gm of Ammonium Molybdate in 100 ml of water

Principle: Inorganic phosphorus reacts with ammonium molybdate in an acidic medium to form a phosphomolybdate complex.



Procedure:

➤ [2-ml of urine] + [0.5 ml concentrated HNO₃] + [3 ml of 5% Ammonium Molybdate], Heat

Canary yellow precipitate of Ammonium phosphomolybdate are formed


Sulphate:

Reagent:

1 % HCL : Take 1 ml of concentrated HCL & make upto 100 ml

10% Barium chloride :Dissolve 50 gm of Barium chloride in 500 ml of water

Principle:



Procedure:

[5 ml urine] + [1 ml 1 % diluted HCL] + [2 ml of 10% Barium chloride].
White precipitate of BaSO₄ are formed

Organic Chemical constituents


Urea:

(Specific Urease Test)


Reagent:

1% phenolphthalein :Dissolve 1 gm of phenolphthalein in 100 ml methanol

Principle:



CO₂ evaporates.



In this reaction the liberation of NH₃ changes the pH to alkaline side, turning phenolphthalein to pink colour.

Procedure:

[2 ml Urine] + [2 drops phenolphthalein]

Add 2% Na₂CO₃ till faint pink color is seen.

Add acetic acid, one drop at a time, with mixing, till faint pink color just disappears.

Add a spatula of Urease powder(Jack Bean Meal Powder), mix.
Pink color develops after few minutes.

Uric acid:

Phosphotungstic acid reduction test:

Reagent:

10% Sodium carbonate: Dissolve 10gm of sodium carbonate in 100ml of water

Phosphotungstic acid Reagent :

Stock : Dissolve 50 gm of sodium tungstate in 400 ml of water & add 40ml of 85% phosphoric acid .Make final volume to 500ml.

Working : Dilute 50ml of stock to 500ml with water.

Principle:

Uric acid is reducing agent in alkaline medium.It reduced phosphotungstic acid into tungsten blue.

Procedure

To 2.5 ml of urine add 0.5 ml of sodium carbonate and 0.5 ml of Phosphotungstic acid reagent working reagent

Creatinine

Reagent:

Refer SOP in dokuwiki document

Creatinine R1 (NaOH)

1. Weigh 24 gm NaOH.
2. Dissolve in approximately 500 ml DI water.
3. Add 20 ml of 30% brij in above mixture.
4. Weigh 2 gm SDS and pour it into approximately 200 ml water in beaker. Heat the solution until SDS dissolve.
5. Add SDS containing solution in main mixture.
6. Make upto 2 liter with DI water.

Creatinine R2 (Picric acid)

1. Dry picric acid between filter paper pieces
2. Weight 9.16 gm dry picric acid
3. Dissolve in approx. 600 ml water
4. Add 20 ml of 30% Brij in above mixture.

5. Remove froth with a clean object of glass or plastic dipped in capryl alcohol
6. Make 2 liter with water

Alkaline picrate reagent:

Mix R1 and R2 1:1 to make working alkaline picrate reagent.

Principle:

Creatinine forms creatinine picrate in alkaline medium which is orange in colour

Procedure

2 ml alkaline picrate solution + 1 drop of urine & mix

Physical Characteristics of Urine.

Physical characteristics of urine	Observation	Interference
Volume		
Appearance		
Colour		
Odour		
Reaction		
Specific gravity		

Inorganic constituents of urine.

Inorganic constituents of Urine	Observation	Interference
Ammonia		
Chloride		
Calcium		
Phosphorus		

Sulphate		
----------	--	--

Organic constituents of urine

Organic constituents of Urine	Observation	Interference
Urea		
Uric acid		
Creatinine		

6. Pathological Urine

Appearance:

Turbid: infection (cells make urine turbid)

Color:

Yellow: Hepatic jaundice & obstructive jaundice (Conjugated bilirubin)

Red: Hematuria, rifampicin therapy

Red on exposure to air: porphyria

Black on exposure to air: alkaptonuria

Odour:

Fruity: diabetic ketoacidosis (acetone)

Mousy smell: Phenylketonuria. (Phenylacetyl glutamine)

Foul smell: Urinary tract infections. (H₂S etc.)

Specific Gravity :

High Specific Gravity

- Diabetes mellitus
- Diarrhoea
- Dehydration

Low Specific Gravity

- Diabetes insipidus.
- Renal failure
- Excessive fluid intake
- Acute tubular necrosis

Abnormal Constituents of Urine

I. Protein:

Reagent :

A. Sample preparation :

10 mg% albumin : Dissolve 100 mg bovin albumin in 1000 ml of water

50 mg% albumin: Dissolve 500 mg bovin albumin in 1000 ml of water

100 mg% albumin: Dissolve 1000 mg bovin albumin in 1000 ml of water

B. 1% Acetic acid:

- 5 ml of acetic acid in 500 ml of water

C. 30% Sulphosalicylic acid:

- Dissolve 150 gm of Sulphosalicylic acid in 500 ml of water

Proteinuria and albuminuria

Proteinuria

Normal Adult	<150 mg /day
Proteinuria	>=150 mg /day
Proteinuria	>3500 mg / day

Albuminuria

Normal Adult	<30 mg /day
Microalbuminuria	30-300 mg /day
Macroalbuminuria	>300 mg /day

Albumin (Filtered but not reabsorbed) and Tamm-Horsfall protein (secreted by renal tubules) are normally present.

Causes of Proteinuria:

- Pre-renal: (overload proteinuria)
 - o (Many non-Albumin proteins)
 - o Multiple myeloma (light chains of immunoglobulins)
 - o Severe hemolysis (Hemoglobin)
 - o Severe muscle injury (Myoglobinuria)
- Renal: *Glomerular diseases* (Mainly albumin, being small)
 - o After streptococcal infection
 - o Diabetes mellitus
 - o Hypertension
 - o Lipoid Nephrosis (Nephrotic range proteinuria)
 - o *Tubular diseases* (decreased reabsorption of proteins)
 - o (Small, normally reabsorbed, proteins like Beta2 microglobulin, Retinol Binding protein)
 - o Tubular necrosis due to Drugs and toxins
- Post Renal: (various blood and cellular proteins)
 - o Bleeding in urinary tract
 - o Infection in urinary tract
 - o Tumor in urinary tract
- Other causes:
 - o Postural: on standing posture.
 - o Exposure to cold, physical activity, fever.
 - o Last weeks of pregnancy

Heat coagulation Test:

Principle :

Proteins have net zero charge at their iso-electric pH (pI). So, at pI, protein molecules have minimum repelling force. Thus proteins are easily precipitated at pI.

When proteins are heated, weak bonds like hydrogen-bonds, salt bonds and van-der-wal forces are broken. Proteins are said to be denatured. Core hydrophobic regions of denatured Albumin can form intermolecular associations and cause precipitation. Thus, in order to precipitate proteins like albumin, two conditions are required.

1. Bring albumin to its pI(5.4) by adding few drops of 1% acetic acid
2. Heat the solution

Procedure

Fill 3/4 th of the test tube with urine sample,Heat the upper part on the flame till either turbidity appears or urine starts boiling.Then add few drops of 1% acetic acid if turbidity develops & note change.

In case of multiple myeloma, light chains of immunoglobulin precipitate between 40-60 degrees centigrade. With further heating turbidity disappears. Turbidity appears again on cooling to 40-60 degree centigrade.

Sulphosalicylic Test:

Principle

Test is based on the precipitation of urine protein by a strong acid, sulfosalicylic acid. Precipitation of protein in the sample seen as increasing turbidity)

Unlike the routine urine protein chemistry dipstick pad, the SSA reaction will

detect globulin and Bence-Jones proteins, in addition to albumin

Method: 3 ml of urine + 0.3 ml of 30% Sulphosalicylic acid, mix.

Turbidity indicates presence of urinary proteins.

Iodinated contrast agents used for evaluation of renal disorders can give the test positive.

False positives:

X-ray contrast media

High concentration of antibiotics, such as penicillin and cephalosporin derivatives.

False negatives:

Highly buffered alkaline urine. (The urine may require acidification to a pH of 7.0 before performing the SSA test.)

Dilute urine

Turbid urine - may mask a positive reaction. Again, best practice is to always used supernatant from a properly spun urine sample.


Dip-Stick Test:

Principle

Testing for protein is based on the phenomenon called the "**Protein Error of Indicators**" (ability of protein to alter the color of some acid-base indicators without altering the pH).

This principle is based on the fact that proteins alter the colour of some pH indicators even though the pH of the media remains constant. This occurs because proteins (and particularly albumin) acquire hydrogen

ions at the expense of the indicator as the protein's amino groups are highly efficient acceptors of H⁺ ions.



At pH 3 and in the absence of proteins both indicators are yellow, as protein concentration increases the colour changes through various shades of green until it becomes a dark blue.

According to the manufacturer, the strip's protein pad contains tetrabromophenol blue or 3',3',5',5'-tetrachlorophenol-3,4,5,5-tetrabromosulphonphthalein, as well as an acid buffer substance to maintain pH at a constant level.

The main problem with the protein tests found on urine test strips is that very alkali urine can neutralise the acid buffer and produce a false positive reading that is unrelated to the presence of proteins. Another similar error occurs if the strip is left submerged in the urine sample for too long.

This method is more sensitive to albumin than to globulin, Bence Jones protein and mucoprotein are examples of globulin components that are sometimes present in urine, but are not distinguishable by the dipstick method for protein

Method: Dip the strip for Albumin in urine. Drain excess urine from strip. Read the color chart. (Read instruction manual provided with the strips for time of reading after dip.).

Because the dipstick test detect albumin, it can not identify many pre-renal proteinuria caused by Hb, Mb and light chains of Igs.

All the three tests mentioned above are qualitative and used for screening proteinuria and albuminuria. Once proteinuria is found quantitative estimation of proteinuria and albuminuria is required for clinical decision making.

What Will You Do:

Perform all three tests with urine. Draw table of your observations.

Sr. no	Concentration	Heat coagulation test	Dipstick test	Sulphosalicylic acid	Interference
1	10 mg %				

2	50 mg %				
3	100 mg %				
4	Urine sample				

Which of the three tests is most sensitive?

Write biochemical explanation of proteinuria in diabetes mellitus and hypertension.

II. Acetone & acetoacetic acid (Ketone Bodies):

Reagent

1. Ammonium sulphate powder
2. Small crystals of sodium nitroprusside
3. liquor Ammonia

Rothera's powdered reagent :

- Sodium Nitroprusside 1 gm
- Sodium carbonate 20 gm
- Ammonium sulphate 20 gm

Mix & grind all in fine particales & stored in air-tight container.

Sample Preparation

0.1 ml/L Acetone : Take 0.1 ml Acetone in 1000 ml DI water

1 ml/L Acetone : Take 1 ml of Acetone in 1000 ml DI water

10 ml/L Acetone :Take 10 ml of Acetone in 1000 ml DI water

Principle

Acetoacetic acid and acetone form a violet coloured complex with sodium nitroprusside in alkaline medium. Acetoacetic acid reacts more sensitively than acetone. Values of 10 mg/dl of acetoacetic acid or 50 mg/dl acetone are indicated. Phenylketones in higher concentrations interfere with the test, and will produce deviating colours. β -hydroxybutyric acid (not a ketone) is not detected.

Sodium Nitroprusside : acetone form a violet coloured complex with sodium nitroprusside in alkaline medium

Sodium carbonate: Provide Alkaline medium

Ammonium sulphate : Precipitate other protein which give purple colour with sodium nitroprusside & make solution Heavier than liquire Ammonia, so Ammonia may be remain on top of solution ,so purple ring is formed.

Rothera's test, liquid reagent

Saturate 2ml urine with ammonium sulphate powder.

Add a small crystal of sodium nitroprusside. Mix.

Add 0.5 ml liquor ammonia by side of the tube to form a ring.

Permanganate/Purple color ring is formed

Rothera's test, powdered reagent

Take a pinch of Rothera's powdered reagent

Add 1-2 drops of urine on powder.

Permanganate/purple color is formed

What Will You Do:

Perform both tests with given sample of urine.

Perform both tests with 0.1ml/L , 1/ml/L , 10 ml/L acetone

Sr. no	Concentration	Rothera's test,powder reagent	Rothera's test, Liquid reagent	Interference
1	0.1 ml/L			
2	1 ml/L			
3	10 ml/L			
4	Urine sample			

Perform both tests with given sample of urine.

Perform both tests with 0.1ml/L , 1/ml/L , 10 ml/L Acetoacetate

Sr. no	Concentration	Rothera's test,powder reagent	Rothera's test, Liquid reagent	Interference
1	0.1 ml/L			
2	1 ml/L			
3	10 ml/L			

4	Urine sample			
---	--------------	--	--	--

Which other tests in blood and urine are usually done when tests for ketone bodies are positive?

III. Bile Salts:

REAGENT

Bile salt sample : Dissolve 2 gm of Bile salt powder into 1000 ml of water.
Sulfur powder

Principle

Sulphur powder is non-polar. It floats on water surface due to surface tension of water. Bile salt reduces surface tension of water and thereby sulphur powder sinks.

Procedure

Hay's sulfur flower Test:

Sprinkle a pinch of sulphur powder over 2 ml urine in a test tube & Sprinkle a pinch of sulphur powder over 2 ml Water in a test tube. Observed & compare immediately without shaking of test tubes. Sulphur powder sink to the bottom of the test tube if bile salts are present.

What Will You Do:

Perform the Hay's sulfur flower test with given sample

Sample	Observation	Interference
Bile salt solution		
Water		
Urine sample		

IV : Glucose

Perform both the tests with urine. Draw table of your observation.

Perform both tests with 100 mg%, 500 mg%, 1 gm% glucose. Note color of the test. Draw table showing the results as follows.

Glucose %	Benedict's Test color	GOD Strip test color	Interference
100 mg%			
500 mg%			
1000 mg%			
Urine sample			

7. Estimation of acid output by stomach.

Parietal cells of gastric mucosa secrete H^+ using $H^+-K^+-ATPase$. Gastrin, acetylcholine (from vagus) and histamine stimulate H^+ secretion. Thus, abnormality of *parietal cells*, *G cells* and *Vegas* are important in disturbances of gastric acid secretion.

Collection of Gastric juice for Analysis

Preparation of the Patient:

1. Patient should be fasting 10 to 12 hours (preferably since bed time night before).
2. Patient should have not received any medications specially anticholinergic agent, H_2 blocker, Antacids since night before as they are liable to alter the results.
3. Procedure should be explained to the patient in simple words.

Procedure:

1. Remove dentures if there are any from patient's mouth. Take the Nasogastric tube and lubricate it.
2. Check patient's nostrils and choose the one (while patient still sitting up right with neck flexed) through which breathing is easier and nostril is wider.
3. Begin intubation by gently pushing the tube, some times tube gets-curved up in the pharynx and there is excessive coughing or gagging which may prevent further passage, at this time tube is drawn back a few inches, patient, is reassured and intubation is resumed.
4. During this insertion patient is instructed to swallow and continue to swallow throughout the intubation period.
5. After the tube has progressed to approximately 40 cm (the first mark on the tube), the head may be allowed to resume its comfortable position.
6. Continue intubation by gently pushing the tube with the patient still swallowing until the fourth mark or 65 cm is reached.
7. Tape the tube to the patient's nose with adhesive tape. At this point patient is sent to the X-ray for fluoroscopy to check position.
8. The tube should lie along the lesser curvature with the tip in the antrum of the stomach.
9. In patient with partial gastrectomy tip of tube should be in the most dependent portion of the stomach.

Collection of Gastric Juice For Analysis

1. Empty the stomach of its contents with a 50 cc. Syringe.
2. After recording the pH volume and colour, this residual volume may be discarded.
3. After emptying stomach of the residual volume, collection of gastric juice is begun under Basal conditions.
4. At least four samples are collected each 15 minutes apart in separate containers.

5. Collection may be carried out either manually with the syringe or by using a suction pump.
6. During this procedure patency of nasogastric tube is maintained by injecting about 50 cc of air down the tube.
7. Gastric fluid specimen should be spot checked as a guide to whether the patient is making-acid or is achlorhydric.
8. After having collected gastric juice under basal conditions augmented or stimulated gastric analysis may be carried out as follows:
9. Pentagastrin administered by sub-cutaneous injection in the dose 6 mg per kg body weight. It is a synthetic peptide having the biologically active sequence of gastrin.
10. The gastric secretion is collected every 15 minutes for next 1 hour.

1. Basal acid output (B.A.O):

It is the acid output in milimol per hour in the basal secretion.

2. Maximal acid output (M.A.O)

It is the acid output in milimol per hour, given by the sum of acid output of the four 15 minute sample after the stimulation

Hypochlorhydria: (decreased acid output, pH>4)

- Pernicious anemia
- Autoimmunity to parietal cells destroys them.
- Antibodies to Na⁺-K⁺-ATPase are found
- Chronic *Helicobacter Pylori* infection of gastric mucosa.
- Treatment with Proton pump inhibitors, H₂-Blocker
- Vagotomy

Hyperchlorhydria: (increased acid output)

- Zollinger-Ellision Syndrome
- G cells tumors in GIT

Reagent

0.1 mol/L NaOH : Dissolve 20 gm of NaOH in 5000 ml of water

1 % phenolphthalein : Dissolve 1 gm of phenolphthalein in 100 ml methanol

Sample preparation

Gastric juice Sample : 0.1mol/L HCL solution

How 0.1 mol/L Hcl will be prepared?

1000ml of HCL solution contain=11.5 mol H⁺

??????? =0.08 mol H⁺

=1000x0.1/11.5

=8.6 ml

So add 17 ml of concentrated HCL & make upto 2 liter with water.

Examples

Example-1:

If you want your **result** will be Gastric Acid Output (mmol/hr) = 5 mmol/hr and You give Fasting Gastric juice output in 1 hour =100 ml/hr then prepare gastric juice sample as follow,

Fasting Gastric juice output =100 ml/hr
BAO = 5 mmol/L

100 ml of fasting gastric juice contain = 5 mmol/L HCL

1000 ml of fasting gastric juice contain = ???

$$= 1000 \times 5$$

$$100$$

$$= 50 \text{ mmol/L HCL}$$

$$= 0.05 \text{ mol/L HCL}$$

Now We use fixed 10 ml of Gastric juice sample & titrate with fixed 0.1 mol/L NaOH

10 ml of 0.05 mol/L HCL =-----ml of 0.1 mol/L NaOH

V1=10 ml of Gastric juice

V2=???? ml of NaOH

N1=0.05 mol/L HCL

N2=0.1 mol/L NaOH

$$V2=10 \times 0.05/0.1$$

$$=5 \text{ ml of } 0.1 \text{ mol/L NaOH}$$

Thus 5 ml of 0.1 mol/L NaOH is required to titrate 10 ml of 0.05 mol/L HCL.

Now ,Check your sample of gastric juice is made proper or not by following formula,

Gastric Acid Output =

[Average Reading R] *[Gastric Juice Output in one hour]

$$100$$

We require 5 ml of NaOH & give 100 ml/hr Gastric output,so our result is

Gastric Acid Output = 5 x 100/100

$$=5 \text{ mmol/hr ,that is our BAO.}$$

Example-2:

If you want your **result** will be Gastric Acid Output (mmol/hr) = 8 mmol/hr and You give Fasting Gastric juice output in 1 hour =80 ml/hr then prepare gastric juice sample as follow,

Fasting Gastric juice output =80 ml/hr
BAO = 8 mmol/L

80 ml of fasting gastric juice contain = 8 mmol/L HCL
1000 ml of fasting gastric juice contain = ???
= 1000 x 8/80
= 100mmol/L HCL
= 0.1 mol/L HCL

Thus Take 8.6 ml of Concentrated HCL solution and make upto 1000ml with water is made to 8.6 mol/L HCL solution.

Now We use fixed 10 ml of Gastric juice sample & titrate with fixed 0.1 mol/L NaOH

10 ml of 0.1 mol/L HCL = _____ ml of 0.1 mol/L NaOH
V1=10 ml of Gastric juice V2=???? ml of NaOH
N1=0.08 mol/L HCL N2=0.1 mol/L NaOH

V2 = 10 x 0.1/0.1
= 10 ml of 0.1 mol/L NaOH

Thus 10 ml of 0.1 mol/L NaOH is required to titrate 10 ml of 0.1 mol/L HCL.

Now , Check your sample of gastric juice is made proper or not by following formula,

Gastric Acid Output =
[Average Reading R]* [Gastric Juice Output in one hour]

100

We require 10 ml of NaOH & give 80 ml/hr Gastric output, so our result is Gastric Acid Output = 10 x 80/100
= 8 mmol/hr ,that is our BAO.

Principle:

Acid output in stomach is measured as mmol/hour. For its measurement, amount of gastric juice output as well as amount of acid in gastric juice needs to be measured.

Amount of Gastric juice output is measured by suction of gastric juice using Ryle's tube inserted in to stomach.

Amount of acid in gastric juice is measured as follows.

Free Acidity:

Due to H^+ (H_3O^+) ions.

Combined Acidity:

Some of the H^+ in gastric juice are bound to other anions like proteins and lactic acids at low pH of Gastric Juice. These represent combined acidity.

(Proteins $^-$).(H^+) , (Lactate $^-$).(H^+)

Free Acidity + Combined Acidity = Total acidity

On addition of alkali, initially free H^+ and later on combined H^+ are neutralized.

When not much H^+ remain in solution (at pH 8.6), Phenolphthalein indicator becomes pink. The requirement of alkali is used to calculate acid output.

Procedure:

First Reading:

- Step 1
 - o Take 10 ml gastric juice in a flask/beaker.
 - o Add 1 drop of phenolphthalein.
 - o (Do not mouth pipette anything)
- Step 2
 - o Fill burette with 0.1 mol/L NaOH up to zero mark.
 - o Perform as follows.
- Step 3
 - o Add 1 ml of NaOH from burette, mix, and watch for pink color.
 - o Repeat above step till pink color develops.
 - o Suppose reading is X_1 ml of NaOH

Second Reading:

- Repeat-step 1 and step-2 of first reading.
- Add [$X_1 - 1$] ml of NaOH from burette mix.
- Add NaOH drop wise till pink color develops.
- Take Reading will be X_2 .

Third Reading:

- Repeat-step 1 and step-2 of first reading.
- Add [$X_2 - 1$] ml of NaOH from burette mix.
- Add NaOH drop wise till pink color develops.
- Take Reading will be X_3 .

Find Average(**R**) of X_2 & X_3 .

Calculation:

Explanation of calculation:

1 mol NaOH \equiv 1 mol HCl

R ml of 0.1 mol/L NaOH	\equiv R	ml of 0.1 mol/L HCl
	\equiv R / 10	ml of 1 mol/L HCl
	\equiv R / (10*1000)	mol/ ml of HCl
	\equiv (R / 10)	mmol/ml of HCl

Thus, 10 ml of Gastric Juice will have (R/10) mmol HCl equivalents.

Thus, 1 ml of Gastric Juice will have (R/100) mmol HCl equivalents.

Gastric acid output = (R/100) * G mmol/hr

G = Gastric Juice Output (ml/hr)

Normal Gastric Juice Output (ml/hr) = 80 ml/hr

Result:

Gastric Acid Output (mmol/hr) =

Reference Ranges:

- Fasting Gastric Juice Output: 20-100 ml /hr
- Basal Acid Output (BAO): Measured in fasting state
 - o Normal 1-6 mmol/hr
 - o ZE Syndrome >15 mmol/hr (M)
 - >10 mmol/hr (F)
- Maximum Acid Output (MAO): Measured after pentagastrin stimulation
 - o Normal 5-40 mmol/hr

In pernicious anemia, both MAO and BAO are almost zero.

Above reference ranges are not universally accepted. Serum gastrin level, pH of gastric juice and other clinical finding e.g megaloblastic anemia are important to establish diagnosis.

What will you do:

Estimate gastric acid output in given sample or gastric juice.

Consider Gastric juice output 80 ml/hr.

Initial reading (X)	Next reading (Y)	Volume of NaOH used (X-Y) ml
Average of (X-Y)		

Gastric Acid Output =
 [Average Reading R] * [Gastric Juice Output in onehour]

 100

Result : Your Gastric acid output is -----

Comment on your result

Q-1 What is Zollinger-Ellision syndrome?

Q-2 What happens to Gastric acid output in the ZE syndrome? Why?

Q-3 Write complications of the ZE syndrome.

Q-4 Write cause of destruction of parietal cells in pernicious anemia.

Q-5 What happens to Gastric acid output in the pernicious anemia?
Why?

Q-6 Which other important products are formed and secreted by parietal cells?

Q-7 Why should destruction of parietal cell lead to anemia?

Q-8 What is difference between gastrin and pentagastrin.

Q-9 Both pernicious anemia and ZE syndrome result in high serum gastrin level. Explain.

Q-10 Explain mechanism of action and use of ranitidine and omeprazole as drugs.

8. Secretion and buffering of acids by kidney.

Reagent

1. 1 % phenolphthalein : Dissolve 0.5 gm of phenolphthalein in 50 ml of Methanol.
2. Neutral formalin (formaldehyde): Take 500ml of formaldehyde & add 0.1ml of phenolphthalein in solution. Then add 0.1 mol/L NaOH till colorless formaldehyde solution become slight pink coloured.
3. 0.1mol/L NaOH : Dissolve 20 gm of NaOH & make upto 5000 ml with Water.

Urine Sample Preparation:

Urine output ml/day = U

Titration acidity mmol/day = A

$$\text{Take } \frac{A}{U} \times 68 \text{ gm of } \text{KH}_2\text{PO}_4 \quad \text{MW of } \text{KH}_2\text{PO}_4 = 68 \text{ gm/L}$$

Ammonia bound acidity mmol/day = B

$$\text{Take } \frac{B}{U} \times 66 \text{ gm of } (\text{NH}_4)_2\text{SO}_4 \quad \text{MW of } (\text{NH}_4)_2\text{SO}_4 = 132 \text{ gm/L}$$

Here two NH_4^+ is released when 1 molecule of $(\text{NH}_4)_2\text{SO}_4$ will be dissociated.

Example

You want to give Titration acidity = 30 mmol HCL /day &

Ammonia bound acidity = 40 mmol HCL /day ,then prepare Urine sample as follow,



$$\begin{aligned} \text{Urine output U} &= 1500 \text{ ml/day} \\ \text{Titration acidity mmol/day A} &= 30 \text{ mmol HCL/day} \\ &= A/U \times 68 \\ &= 30/1500 \times 68 \\ &= 1.36 \text{ gm of } \text{KH}_2\text{PO}_4 \\ \\ \text{Ammonia Bound acidity mmol/day B} &= 40 \text{ mmol HCL/day} \\ &= B/U \times 66 \\ &= 40/1500 \times 66 \end{aligned}$$

=1.76 gm of $(\text{NH}_4)_2\text{SO}_4$

Finally dissolve 1.36 gm of KH_2PO_4 and 1.76 gm of $(\text{NH}_4)_2\text{SO}_4$ & make upto 1000 ml with water.

Principle:

Catabolism of food substances produces H^+ and OH^- . In the process, there is excess of H^+ over OH^- . Excess H^+ is excreted by kidney. NH_3 and Phosphate buffer the H^+ secreted by renal tubules.




You will estimate total Acids in urine and proportions buffered by ammonia and phosphate.

Correlate the experiment with theoretical concepts of renal regulation of pH learnt in the classroom.

pK of reaction (1) is 9.25.

pK of reaction (2) is 6.8.

For phenolphthalein: pH 8.2 = colorless; pH 10 = red



Phenolphthalein is a weak acid, which can lose H^+ ions in solution. The phenolphthalein molecule (HIn) is colorless, and the phenolphthalein ion (In^-) is pink. When a base is added to the phenolphthalein, the molecule \rightleftharpoons ions equilibrium shifts to the right, leading to more ionization as H^+ ions are removed

When urine, acidic in nature, is titrated with NaOH, initially reaction (2) goes towards left. When all H_2PO_4^- is converted into HPO_4^{2-} , pH rises to 8.6, causing ionization of phenolphthalein. Phenolphthalein ion produced pink colour, so solution turns into pink coloured. NaOH required to reach this stage represents H^+ bound to phosphate, called **“Titrable Acidity”**.

Neutral formalin is added to urine. We will convert formalin (Acid) to Neutral formalin, otherwise formalin (acid) itself reacts with NaOH when we measure H^+ of NH_4^+

Now, Formaldehyde is added to urine. Following reaction occurs.



Released H^+ decrease pH of urine, making phenolphthalein colorless again.

Further titration with NaOH, till phenolphthalein become pink, will actually represent H^+ bound with ammonia released during reaction (3).

It is called "**Ammonia bound acidity**".

H^+ bound to NH_3 can not be titrated without adding formaldehyde.

Hence, H^+ bound to phosphate is called titrable acidity.

Procedure:**First Reading:**

- Step - 1
 - Take 25 ml urine in a flask/beaker.
 - Add 1 drop of phenolphthalein. (Do not mouth pipette anything)
- Step - 2
 - Fill burette with 0.1 mol/L NaOH up to zero mark.
- Step - 3
 - Perform as follows.
 - Add 1 ml of NaOH from burette, mix, and watch for pink color.
 - Repeat above step (adding 1 ml NaOH) till pink color develops.
 - Suppose reading is X_1 ml of NaOH
- Step - 4
 - Add 10 ml of neutral formalin & Mix.
 - The pink color disappears.
 - Repeat step-3.
 - Suppose the reading is Y_1

Second Reading:

- Repeat-step 1 and step-2 of above.
- Add [$X_1 - 1$] ml of NaOH from burette & mix.
- Add NaOH drop wise till pink color develops.
- Take reading X_2 .
- Add 10 ml of neutral formalin & Mix.
- Add [$Y_1 - 1$] ml of NaOH from burette & Mix.
- Add NaOH one drop wise till pink color develops. Take reading Y_2 .

Third Reading:

- Repeat-step 1 and step-2 of above.
- Add [$X_2 - 1$] ml of NaOH from burette & mix.
- Add NaOH drop wise till pink color develops.
- Take reading X_3 .
- Add 10 ml of neutral formalin & Mix.
- Add [$Y_2 - 1$] ml of NaOH from burette & Mix.
- Add NaOH one drop wise till pink color develops. Take reading Y_3 .

Find X (Average of X_2 & X_3) and Y (Average of Y_2 & Y_3) .

Explanation of calculation:**Titration acidity: reading X ml**

1 mol NaOH \equiv 1 mol HCl

$$\begin{aligned}
 X \text{ ml of } 0.1 \text{ mol/L NaOH} & \equiv X \text{ ml of } 0.1 \text{ mol/L HCl} \\
 & \equiv X / 10 \text{ ml of } 1 \text{ mol/L HCl} \\
 & \equiv X / (10 \times 1000) \text{ mol/ml of HCl} \\
 & \equiv (X / 10) \text{ mmol/ml of HCl}
 \end{aligned}$$

As titration is done with 25 ml of urine,

$$\text{Titration acidity in 25 ml of urine} = (X / 10) \text{ mmol HCl}$$

$$\text{Titration acidity in 1 ml of urine} = X / (10 \times 25) \text{ mmol HCl}$$

If urine output per day is U ml

$$\text{Excreted Titration acidity / day} = (U \times X) / 250 \text{ mmol HCl}$$

Ammonia bound acidity: reading Y ml

It is expressed either as **mmol of HCl** or **mg of ammonia**

$$\text{Ammonia bound acidity / day} = (U \times Y) / 250 \text{ mmol HCl}$$



$$1 \text{ mmol of NH}_3 \text{ binds } 1 \text{ mmol of H}^+ \text{ to form } 1 \text{ mmol of NH}_4^+ \quad \text{---(b)}$$

$$\text{MW of Ammonia (NH}_3) = 17 \text{ gm}$$

$$1 \text{ mmol NH}_3 = 17 \text{ mg of NH}_3 \quad \text{---(a)}$$

From (a) and (b)

$$\begin{aligned}
 \text{Excreted Ammonia / day} & = (U \times Y) / 250 \text{ mmol NH}_3 \\
 & = ((U \times Y) / 250) \times (17) \text{ mg NH}_3
 \end{aligned}$$

$$\text{Excreted Ammonia / day} = U \times Y \times (0.068) \text{ mg NH}_3$$

Reference Range:

$$\text{Titration acidity} = 20\text{-}50 \text{ mmol HCl / day}$$

$$\text{Ammonia bound acidity} = [30\text{-}50 \text{ mmol HCl/day}] \text{ or } [510\text{-}850 \text{ mg NH}_3/\text{day}]$$

$$\text{Total acid excretion} = 70\text{-}100 \text{ mmol/day}$$

What will you do:

Estimate Titration and ammonia bound acidity in given sample of urine.

Titration acidity

No.	Initial reading(ml)	Final reading(ml)	Difference(ml)
X₂			

X₃			
Average X			

Ammonia bound acidity

No.	Initial reading(ml)	Final reading(ml)	Difference(ml)
Y₂			
Y₃			
Average Y			

Result & conclusion

Titration acidity =

Ammonia bound acidity =


What is the source of phosphate in urine?

What is the source of ammonia in urine?

Diabetic ketoacidosis elevate urinary ammonia. Explain.

9. Colorimetry

Colored molecule absorbs various wavelength of light passing through their solution.



For a given wavelength of light, ratio of (emerging light intensity) to (imparted light intensity) is called Transmittance T.

$$T = e^{-kct}$$

c = concentration of colored molecule

t = length of light path

k = constant

$$-kct = \log_e T$$

$$-kct = \frac{\log_{10} T}{\log_{10} e}$$

$$-k'ct = \log T \quad (\text{common logarithm}) \quad k' = \text{constant}$$


$$-k'ct = \log (T \cdot 100 / 100)$$

$$-k'ct = \log (T \cdot 100) - 2$$

$$k'ct = 2 - \log (T\%) \quad (T\% \text{ is called percentage Transmittance})$$

$$k'ct = A \quad (2 - \log T\%) \text{ is called Absorbance, denoted as } A$$

Following graph describe relationship between T% and A.



$A = k'ct.$ (Beer's and Lambert's law)

Hence,

$A \propto t$ Absorbance is proportional to length of light path

$A \propto c$ Absorbance is proportional to concentration of substance

Therefore, If light path is constant, for concentration (C_1 and C_2) and respective absorbance (A_1 and A_2)

$$C_1/C_2 = A_1/A_2$$


If A_1 and A_2 is measured and C_2 is known

$C_1 = (A_1/A_2) * C_2$ can be calculated. -----(1)

This principle is utilized by biochemistry laboratory to measure various substances in biological materials.

Various instruments based on the principle are colorimeter and spectrophotometer.

Instrument:



Light source emit light of all wavelengths

Monochromator allow only certain wavelength of light to pass. (Mono + color)

Cuvette is a transparent vessel holding colored solution

I_o = Imparted light

I_t = Immerging light

Photocell converts light in to current. Current is proportional to light intensity.


Galvanometer measures current.

General procedure to use colorimeter:

Suppose concentration of Glucose in plasma is to be estimated.

Glucose is colorless, hence can not be measured directly.

Add fixed amount of Y in fixed amount of plasma. P and Q are produced



Suppose Q is colored compound and absorbs light of a particular wavelength. Its concentration will be proportional to concentration of glucose.

Take a solution of glucose with known concentration C (it is called calibrator) and process as above in 2.

Take a water (it is called blank) and process as above in 2.

Measure absorbance of color produced by Sample and Standard and blank.

Blank Absorbance, amount of color produced with no glucose, needs to be deducted from absorbance of Sample and Standard.

Using equation (1)

$$\text{Glucose concentration in plasma (mg\%)} = \frac{(A_{\text{sample}} - A_{\text{blank}})}{(A_{\text{standard}} - A_{\text{blank}})} * C$$

What will you do:

Reagent:

Buffer:

pH	Chemical drug	Mol/L	MW	Gm/L
6.853	Na2HPO4	0.025	141.96	3.549
	KH2PO4	0.025	136.09	3.402
9.139	Na2 tetraborate	0.01	381.37	3.814

Stock 1

Red coloured solution: Dissolve 20 mg of Phenol red in 50 ml of 9.139 pH buffer

Stock 2

Blue coloured solution: Dissolve 20 mg of BCG(Bromocresol green) in 50 ml of 6.8 pH buffer.

Note: Dilution of stock coloured solution will be always done with respective Buffer.

Diluted Red Dye from Stock-1:

Dilute 1:30 times of red coloured stock-1 solution with buffer of pH=9.1 by Adding 20 ml of stock-1 red coloured solution into 580 ml of Buffer of pH=9.1

Diluted Blue Dye from Stock-2:

Dilute 1:20 times of Blue coloured stock solution with buffer of pH=6.8 by Adding 10 ml of stock Blue coloured solution into 290 ml of Buffer of pH=6.8

Exercise:1

You will be given a concentrated colored solution.

Dilute it in a series of test tubes as follows. Measure absorbance.

Test tube	Diluted Red Dye from Stock-1 (micro-liter)	pH=9.139 Buffer(pH=9.139) (micro-liter)	Absorbance (A) on 505 nm Filter
0	0	1000	
1	200	800	
2	400	600	
3	600	400	
4	800	200	
5	1000	0	

Draw Graph of various Dilution of dye versus its absorbance

Result & Conclusion:

Exerscise-2**[1]: Diluted Red Dye from Stock-1: (Phenol red dye)**

Measure Absorbance of this red coloured solution on different filters.

Filters(nm)	Absorbance
340	
405	
450	
505	
546	
578	
630	
670	

Absorbance spectrum of phenol red dye solution (Red coloured)

[2]: Diluted Blue Dye from Stock-2: [Bromocresol green dye]

Measure Absorbance of these diluted Blue coloured solution on different filters.

Filters(nm)	Absorbance
-------------	------------

340	
405	
450	
505	
546	
578	
630	
670	

Absorbance spectrum of BCG(Bromocresol green) dye solution (Blue coloured)


Draw Graph of various filter versus absorption on that filter for red coloured solution & Blue coloured solution.

Note : Different colored solutions absorb light at different wavelengths in different proportions.

10. Estimation of serum creatinine

Creatinine is produced from creatine present mainly in muscles. It is filtered by glomerulus of kidney.

Principle:



Red colored Creatinine-Picrate complex, also called Janovaski complex, is measured at 505 nm.

The rate of reaction is proportional to concentration of creatinine.

The rate of reaction is also indicated by rate of rise in Absorbance (ΔA)

Thus, [creatinine] \propto ΔA

ΔA is difference of absorbances in 60 seconds of reaction

ΔA for standard is $\Delta A_{\text{standard}}$ and ΔA for sample is ΔA_{sample}

$$\text{Conc. of Creatinine in sample} = \frac{\Delta A_{\text{sample}}}{\Delta A_{\text{standard}}} \times \text{Standard}$$

Reagents

The timed measurements of Absorbance require sophisticated colorimeters with flow-through cuvette. The reaction mixture is aspirated in the cuvette and Absorbance is measured at different time.

The Laboratory technologist will help to carry out following steps:

- NaOH solution :Refer to SOP for creatinine reagent
- Picric acid solution : Refer to SOP for creatinine reagent

Creatinine Standard

2 mg/dl Creatinine :

Dissolve 0.010 gm of creatinine powder in 500 ml of 0.1 mol/L HCl solution

Creatinine Test sample.

4 mg/dl Creatinine :

Dissolve 0.020 gm of creatinine powder in 500 ml of 0.1 mol/L HCl solution

6 mg/dl creatinine :

Dissolve 0.030 gm of creatinine powder in 500 ml of 0.1 mol/L HCl solution

0.1mol/l HCL solution

Add 17.4 ml of Conc. HCL solution (11.5 molar Conc.HCL solu.) & make

upto 2000 ml with DI water.

Creatinine R 1(NAOH)

1. Weigh 12gmNaOH.
2. Dissolve in approximately 500 ml DI water.
3. Add 10ml of 30% brij in above mixture.
4. Weigh 1gm SDS and pour it into approximately 100 ml water in beaker.Heat the solution until SDS dissolve.
5. Add SDS containing solution in main mixture.
6. Make upto 1liter with DI water.

Creatinine R2 (picric acid)

Dry picric acid between filter paper pieces

1. Weight 4.58gm dry picric acid.
2. Dissolve in approx. 300 ml water
3. Add 10 ml of 30% Brij in above mixture.
4. Remove froth with a clean object of glass or plastic dipped in capryl alcohol
5. Make 1 liter with water.

Working alkaline -picrate reagent:

Mix 50 ml R1 & 50 ml R2 on the day of practical for 50 student.

Procedure

For sample and standard perform following.

1 ml of Alkaline Picrate Reagent

+

0.1 ml sample.

Mix it and Analyzed sample immediately at 505 nm wave length.

Note the Absorption (Optical Density=O.D.) at starting of reaction and at end of reaction (after 60 second)

Calculate : Change in O.D. in 60 second = $\Delta A = A_0 - A_{60}$

Calculation and Result:

$$\text{Conc. of Creatinine in sample} = \frac{\Delta A_{\text{sample}}}{\Delta A_{\text{standard}}} \times \text{Standard}$$

Your result will be -----

Comment on your result :

Reference Range:

Male 0.7 - 1.3 mg%

Female 0.6 - 1.2 mg%

1 mmol = 1000 micromole

Creatinine Molecular weight = 113.12

What will you do:


Express adult plasma creatinine reference range in **micromole/L**.

Write clinical conditions affecting plasma creatinine concentration.

11. Estimation of plasma glucose

Reagent:

Glucose reagent : Dissolve 100mg of 4- Aminoantipyrine dye in 1000ml of DI water and add 1 ml of phenol saturated water .



Note :

Wear goggles & Glove while taking phenol.

Senior person must be present.



Glucose test sample :

Add 3 ml of analytical grade Sodium Hypochlorite solution & make upto 10 ml with DI water

Glucose standard sample :

Add 1 ml of analytical grade Sodium Hypochlorite solution & make upto 10 ml with DI water.

Principle:



Procedure

Reagents	Blank	standard	Plasma
H ₂ O	0.1 ml		
Glucose Standard		0.1 ml	
Plasma			0.1 ml
Glucose oxidase + Peroxidase Reagent (GOD POD reagent)	1 ml	1 ml	1 ml
Mix & Incubate at room temperature for 30 min. Read absorbance at 505 nm			
Absorbance	A _{blank}	A _{standard}	A _{plasma}

Calculation:

$$\text{Glucose concentration in plasma} = \frac{(A_{\text{plasma}} - A_{\text{blank}})}{(A_{\text{std.}} - A_{\text{blank}})} * \text{Standard}$$

A_{blank} :

A_{standard} :

A_{plasma} :

Standard conc.:-----

Glucose concentration in plasma =

Your result will be

Comment on your result:

Reference Ranges:

Fasting Plasma Glucose	Interpretation	Oral Glucose Tolerance	Interpretation
<=110 mg%	Normal	<139 mg%	Normal
111-125 mg%	Impaired Fasting Glucose	140-199 mg%	Impaired Glucose Tolerance
>=126 mg%	Diabetes mellitus	>=200 mg%	Diabetes mellitus

Fasting = No food intake for at least 8 hours

Oral Glucose Tolerance = 75 gm glucose orally after 8 hrs of fasting.

Above results are valid if found on two or more occasions.

What will you do:

Measure Glucose concentration in given sample of plasma.

Draw above table again with **mmol/L** format. (Glucose MW=180 gm).

12. Estimation of serum cholesterol

Reagent

Cholesterol reagent :

Dissolve 100mg of 4- Aminophenabenzene & 1 ml of phenol saturated water and make upto 1000ml with DI water.

Cholesterol test sample :

Add 3 ml of HOCL(analytical grade) & make upto 10 ml with DI water

Cholesterol standard sample :

Add 1 ml of HOCL(analytical grade) & make upto 10 ml with DI water

Principle:

Cholesterol ester $\xrightarrow{\text{Cholesterol esterase}}$ Cholesterol + Fatty acid

Cholesterol + O₂ $\xrightarrow{\text{Cholesterol Oxidase}}$ Cholest-4-en-3-one + H₂O₂

H₂O₂ + 4-aminophenazone + phenol $\xrightarrow{\text{Peroxidase}}$ Quinonamine

Reagents and procedure:

Reagents	Blank	Standard	Plasma
H ₂ O	0.1 ml		
Cholesterol Standard		0.1 ml	
Plasma			0.1 ml
Cholesterol oxidase + Peroxidase Reagent (COD POD reagent)	1 ml	1 ml	1 ml
Mix, incubate at room temperature for 30 min. Read absorbance at 505 nm			
Absorbance	A _{blank}	A _{standard}	A _{plasma}

Calculation :

$$\text{Cholesterol concentration in plasma} = \frac{(A_{\text{plasma}} - A_{\text{blank}})}{(A_{\text{standard}} - A_{\text{blank}})} * \text{Standard}$$

Result:

A_{blank} : _____

A_{std} : _____

A_{plasma} : _____

Standard Conc. _____

Glucose concentration in plasma =

Your result will be -----

Comment on your result:

Reference Ranges:

Desirable: <200 mg/dL

Borderline: 200-239 mg/dL

High: ≥ 240 mg/dL

What will you do:

Measure Cholesterol concentration in given sample of plasma.

Rewrite reference ranges in **mmol/L** format. Cholesterol MW = 386.64 gm

13. Estimation of Serum Total bilirubin

Reagent:

Refer to SOP for bilirubine reagent:

R1 (Caffeine)

1. Dissolved 75 gm caffeine in 900 ml deionised water with constant mixing
2. Add 112 gm Na Benzoate in above mixer with constant mixing
3. Add 112 gm anhydrous Na Acetate in above mixer with constant mixing
4. Add 2 gm disodium EDTA in above mixer with constant mixing
5. Make upto **2 liter** with deionised water
6. Filter if turbid
7. Store in glass container in freeze
8. If crystalline precipitation are seen at 2-8'C, bring solution to room temperature to redissolve it before use

Diazo A

1. Dissolve 10 gm sulfanilic acid in 900 ml deionised water
2. Add 30 ml concentrated HCL in above mixer
3. Make upto **2 liter** with deionised water
4. Store in glass container

Diazo B

1. Dissolve 1.25 gm Na nitrite(NaNO_2) and **make upto 250 ml** with deionised water.
2. Store in brown glass container

R2 (Diazo Mix)

Mix Diazo A 10 ml & Diazo B 0.3 ml

Billirubin test solution


Dissolve 2 mg of billirubin powder & make upto 100 ml with DI water

Billirubin test solution :

Dissolve 4 mg of billirubin powder & make upto 100 ml with DI water

Principle

One molecule of bilirubin reacts with two molecules of diazotized sulfanilic acid (dialo mix) in an acid solution to form two purple azobilirubin molecules (560 nm). Direct bilirubin reacts in water as well as with accelerator (e.g. caffeine, methanol), while indirect bilirubin react only in presence of accelerator.



Reagents

Reagents	Test Blank	Test	Standard Blank	Standard
Sample	0.1 ml	0.1 ml		
Standard			0.1 ml	0.1 ml
R1 (Caffeine Reagent)	1 ml	1 ml	1 ml	1 ml
Incubate for 10 minute				
R2 (Diazo Mix)		0.2 ml		0.2 ml
Diazo Blank Reagent (Diazo A)	0.2 ml		0.2 ml	

Absorbance	A _{Test Blank}	A _{Test}	A _{Standard Blank}	A _{Standard}
------------	-------------------------	-------------------	-----------------------------	-----------------------

Blanks are taken to subtract absorbance caused by hemolysis (resulting in presence of red color of hemoglobin in serum).

Diazo blank reagent does not have sodium nitrite, hence do not produce azobilirubin.

Calculation:

$$\text{Total Bilirubin (mg/dL)} = \frac{(A_{\text{Test}} - A_{\text{Test Blank}})}{(A_{\text{Standard}} - A_{\text{Standard Blank}})} \times \text{Std Conc.}$$

Result

Your result will be-----

Comment

Reference ranges: (For Adults)

- Total Bilirubin 0.2-1.2 mg/dl
- Direct Bilirubin 0.1-0.4 mg/dl
- Indirect Bilirubin 0.2-0.7 mg/dl
- Bilirubin = MW 584.67 gm
- 1 mmol=1000 micromole

What will you do:

Enumerate causes of unconjugated hyperbilirubinemia and mixed hyperbilirubinemia.

Express Reference ranges in micromole/Liter.

Sample for bilirubin should not be exposed to light. Phototherapy is used in treatment of neonatal jaundice. Explain and correlate.

14. Estimation of serum total protein

Except **immunoglobulins**, majorities of plasma proteins are synthesized by **liver**. Various tissues **catabolize** plasma proteins. Plasma protein concentration reflects balance between their synthesis and catabolism. Under certain conditions intact proteins from plasma are also lost through **GIT, urine** and **skin**. Proteins from intravascular compartment may reach other body compartments. Protein concentration may also be affected by change in plasma water.

Reagent:

Refer to SOP for total protein.

1. Weight 3 gm $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$.
2. Dissolve in approx. 500 ml water.
3. Weight 9 gm (Na K Tartrate). $(4\text{H}_2\text{O})$ and 5 gm KI.
4. Add sequentially 9 gm (Na K Tartrate). $(4\text{H}_2\text{O})$ and 5 gm KI in copper sulphate solution.
5. Weight 24 gm NaOH.
6. Add slowly with mixing 24 gm NaOH in 100ml of water.
7. Add slowly with mixing NaOH solution in copper sulphate solution.
8. Make upto 1 liter with water

Total Protein Standard:

Dilute Serum pool with DI water(1:20 ratio)

Prepare 10 ml of standard (0.5 ml pool serum+ 9.5 ml DI water)

Total Protein Test:


Dilute Serum pool with Nine part of DI water(1:10 ratio)

Prepare 10 ml of Test (1 ml pool serum+ 9 ml DI water)

Principle :

Two or more peptide bonds of proteins form coordination complex with one Cu^{2+} in alkaline solutions to form a colored product. The absorbance of the product is determined spectrophotometrically at 540 nm.

Cu²⁺ - peptide complex



Procedure:

Reagents	Blank	Standard	Sample
H ₂ O	0.1 ml	-	-
Protein standard	-	0.1 ml	-
Sample	-	-	0.1 ml
Biuret reagent	1 ml	1 ml	1 ml
Mix and incubate at 37° C temperature for 30 min. Read Absorbance at 540 nm			
Absorbance	A _{blank}	A _{standard.}	A _{sample}

Calculation and result:

$$\text{Total Protein concentration in plasma} = \frac{(A_{\text{serum}} - A_{\text{blank}})}{(A_{\text{standard}} - A_{\text{blank}})} * \text{Standard}$$

Your result will be -----

Comment :

Reference ranges:	Serum proteins	6.0-8.0 g/dL
	Albumin	3.5-5.5 g/dL
	Globulins	2.0-3.6 g/dL
	Fibrinogen	0.2-0.6 g/dL

What will you do:

Q-1 Serum protein reference ranges are lower than that of plasma.
Explain.

Q-2 Why reference ranges for plasma proteins can not be expressed in mmol?

QA-3 Enumerate conditions affecting plasma protein level.

15. Estimation of serum albumin

Different disorders affect different plasma proteins differently. Thus, it is useful to know albumin and globulin concentration in serum, in addition to total protein. Once total protein and albumin (as shown below) are estimated, serum globulin can be calculated.

Reagent :

Albumin Standard: Dilute Serum pool with DI water (1:20 ratio). Prepare 10 ml of standard (0.5 ml pool serum + 9.5 ml DI water).

Total Protein Test: Dilute Serum pool with Nine part of DI water (1:10 ratio). Prepare 10 ml of Test (1 ml pool serum + 9 ml DI water).

BCG reagent : Refer SOP for Albumin reagent preparation.

Add 42mg BCG (MW=698) in approx. 250 ml DI water.

Add 5.9 gm succinic acid (MW=118.09, pK_{A1}=4.2, pK_{A2} = 5.6) in above mixer while constantly mixing.

Add 1.8 ml of 30% Brij-35 In above mixer while constantly mixing.


Add 1 gm NaOH in above mixer while constantly mixing

Add 200 mg Na azide in above mixer while constantly mixing.

If required adjust pH to 4.2

Make upto 1000 ml with volumetric flask with deionised water.

Cu²⁺ - peptide complex



Principle:

BCG = BromoCresol Green



At pH 4.2 BCG is yellowish, while Albumin⁺ BCG⁻ complex is greenish. The green color is measured at 630 nm.

Procedure:

Reagents	Blank	Standard	Sample
H ₂ O	0.1 ml	-	-
Standard	-	0.1 ml	-
Sample	-	-	0.1 ml
BCG reagent	1 ml	1 ml	1 ml
Mix, and read immediately at 630 nm.			
Absorbance	A _{blank}	A _{standard}	A _{sample}

Calculation and result:

$$\text{Albumin concentration in serum} = \frac{(A_{\text{sample}} - A_{\text{blank}})}{(A_{\text{standard}} - A_{\text{blank}})} * \text{Standard}$$

Your Result will be -----

Comment on your result

Reference ranges:

Serum Total proteins	6.0-8.0 g/dL
Serum Albumin	3.5-5.5 g/dL
Serum Globulins	2.0-3.6 g/dL
Plasma Fibrinogen	0.2-0.6 g/dL

What will you do:

Calculate S.Globulin from Sample's Albumin conc. & If Total Protein = 7 gm/dl.

Enumerate conditions of alter A:G Ratio.

16. Estimation of Cerebrospinal fluid protein.

Cerebrospinal fluid is not freely permeable to plasma proteins. Hence, its concentration is almost 1/100 times the plasma. Some proteins are synthesized by the pia matter itself. Under various **CNS inflammatory conditions**, CSF protein is increased due to increased permeability of pia matter as well as due to increased synthesis by it.

Reagent:

Pyrogallol Red reagent: Refer to SOP for Pyrogallol red reagent

PR(Pyrogallol red)

Making Reagent

1. Dissolve pyrogallol red 60 mg in 100 ml of methanol.
2. Store in plastic container.

MB(molybdate)

Making Reagent

1. Dissolve disodium molybdate 0.24 gm in 100 ml of deionized water.
2. Store in plastic container.

Final microprotein Reagent

Making Reagent

1. Dissolve succinic acid 5.9 gm in 900 ml of deionized water.
2. Add sodium oxalate 0.14 gm in above mixture with constantly mixing.
3. Add sodium benzoate 0.5 gm in above mixture with constantly mixing.
4. Add PR(Pyrogallol red) 40 ml in above mixture with constantly mixing. Discard other 60 ml PR(Pyrogallol red).
5. Add (molybdate) 4 ml in above mixture with constantly mixing. Discard other 96 ml (molybdate).
6. Calibrate PH meter and if required adjust PH to 2.5.
7. Make up to above mixture 1 L with deionized water.

CSF Protein Calibrator: Take 0.02 ml of serum protein & make up to 10 ml with DI water

CSF protein Sample : Take 0.04 ml of serum Protein & make up to 10 ml with DI Water

Principle:

pyrogallol red-molybdate complex combine with protein and give colour which is measure at 630 nm.

Procedure:

Reagents	Blank	Standard	Sample
H ₂ O	0.1 ml	-	-
CSF Protein Standard	-	0.1 ml	-
CSF Sample	-	-	0.1 ml
Pyrogallol Red reagent	1 ml	1 ml	1 ml
Mix, wait for 10 min, mix before reading at 630 nm.			
Absorbance	A _{blank}	A _{standard}	A _{sample}

$$(A_{\text{sample}} - A_{\text{blank}})$$

Protein concentration in CSF = $\frac{\text{-----}}{(A_{\text{standard}} - A_{\text{blank}})} * \text{Standard}$

Your result & comment

Reference ranges: 15-45 mg%

What will you do:

Enumerate conditions affecting CSF protein level.

17. Estimation of plasma uric acid

Uric acid is formed by catabolism of **purines**. Uric acid is excreted by **kidney**.

Reagent:

Uric acid test sample :



Add 1.5 ml of analytical grade Sodium Hypochlorite solution & make up to 10 ml with DI water.

Uric acid standard sample :

Add 0.5 ml of analytical grade Sodium Hypochlorite solution & make up to 10 ml with DI water.

Principle:

Uric acid yields allantoin and H₂O₂ on action by uricase. Peroxidase use hydrogen peroxide to oxidize various colorless dyes to red colored quinonimine like dyes measured at 505 nm by absorption photometry .



Procedure:

	Test	Standard	Blank
Serum	0.1 ml		
Standard		0.1 ml	
Water			0.1 ml
Reagent	1 ml	1 ml	1 ml
Measure absorbance at 505 nm			
Absorbance	A _{Serum}	A _{Std.}	A _{Blank}

Calculation

$$\text{Uric acid concentration in Sample} = \frac{(A_{\text{Serum}} - A_{\text{blank}})}{(A_{\text{standar.}} - A_{\text{blank}})} * \text{Standard}$$

Result & Comment:

Reference ranges:

Male : 3.6 - 7.7 mg/dL (214 to 458 micromole/L)

Female : 2.5 - 6.8 mg/dL (149 to 405 micromole/L)

What will you do:

Enumerate conditions affecting plasma uric acid level.

Calculate molecular weight of uric acid from reference ranges given.

18. Electrophoresis

Reagent: Refer to Sop for Serum & Hb electrophoresis

Principle:

Electrophoresis is a refers to the migration of charged molecules under electrical field.

Procedure:

Prepare thin 1 % Agarose gel in appropriate buffer.

Apply appropriate sample in thin line over agarose gel.

Keep gel with sample applied in electrophoretic chamber & connect the gel with buffer through strips of filter paper and apply appropriate voltage.

After sample run is completed, switch off the power supply and remove slide from chamber.

Denature proteins in methanol and dry the gel with heating.

Stain slide with appropriate stain.

Buffer :

Barbiturate Buffer (for protein electrophoresis)

Tris Buffer (for protein electrophoresis)

TEB Buffer (for Haemoglobin electrophoresis)


Supporting Media:

Whatman Filter Paper

Cellulose Acetate Paper

Agarose Gel

Polyacrylamide Gel



Clinical Applications:

Diagnosis of sickle trait and sickle disease.

Diagnosis of multiple myeloma

Questions:

What is agarose? Why it is used to prepare gel.

Name other electrophoresis support media.

How much agarose is required to prepare 15 ml of 1% agarose?

Which sample is used for hemoglobin electrophoresis?

What was voltage, current and duration of electrophoresis demonstrated to you?

What are major hazard of electrophoresis procedure?

What precautions must be taken to avoid them?

What is difference between electrophorogram of serum protein and plasma protein?

Name stain used during demonstration.

Draw hemoglobin electrophoretic patten in normal, sickle trait and sickle cell disease patients. Explain its biochemical basis.

Draw hemoglobin electrophoretic pattern in HbC and HbD carrier patients. Explain its biochemical basis.

Draw serum protein electrophoretic pattern in multiple myeloma. Explain its biochemical basis.

What is monoclonal antibody?

19. Chromatography

Principle : Chromatography is a process in which components of a mixture are separated by differential distribution between a mobile phase and a stationary phase. Components with greater distribution into the stationary phase are retained and move through the system more slowly.

Requirements :

Amino acid standard : 1% amino acid standard

Mobile phase : 12(Butanol):3(Glacial acetic acid):5 (Di water)

Sample Type : Serum,Urine .

Equipments & consumables :Chromatography chamber(air tight), Glass rod, clips, Whatman filter paper, Gloves, pencil, scale, centrifuge, pipettes

Stain : Ninhydrin solution (0.25 %)(250 mg of Ninhydrin powder in 100 ml of methanol/acetone)

Procedure :

Clean hands thoroughly with soap.

Wear gloves before handling filter paper.

Take a Whatman filter paper ,make a horizontal line at one end of filter paper,around 1.5-2 cm above from the edge of the paper.

Put marking at 3.5 cm apart for each sample for sample application.

Repeat sample & standard application for twice once previous sample gets dried .

Take 500 ml of mobile phase reagent in reagent chamber.

Clip the dried filter paper on glass rod,Make sure that distance between each sample & rod is equal.

Put glass rod in chromatography chamber, make sure that sample application sites do not get dipped in the solvent.

Close the chamber air tight .Note the time & allow the separation for 4 hours.

Remove the paper from chamber after 4 hours, allow the paper to dry at room temperature.

Take 0.25% Ninhydrin solution in shallow plastic container big enough to accommodate the entire filter paper. Dip paper in it for few seconds.

Put the paper in incubator at 100°C for 20-25 minute/ till purple bands are seen

Preserve the paper in dark room for latter use.

Calculate Rf value

$$R_f = \frac{\text{Distance from application point to solute center}}{\text{Distance from application point to solvent front}}$$

QUESTIONS :

Name stationary phase in the experiment. Is it mainly hydrophobic or mainly hydrophilic? Explain.

Name mobile phase in the experiment. Is it mainly hydrophobic or mainly hydrophilic? Explain.

List hydrophobic amino acids used in the experiment.

List hydrophilic amino acids used in the experiment.

Comment why some amino acids move faster and other slower during the chromatography.

Why wearing gloves is essential in the experiment?

Why wearing protective eye glass is essential in the experiment?

Name few conditions where abnormal amount of some amino acids are lost in urine. Explain their biochemical basis.

Clinical Case - 1

Early in the morning, 40 years old male patient came in emergency with complain of chest pain, perspiration and altered consciousness for 4 hours. Patient also had diabetes mellitus for 10 years. He was taking medicine for diabetes mellitus irregularly. In history, it was found that he was chronic alcoholic and a day before chest pain, he also had heavy alcohol ingestion., with no feed intake. Doctor asked for few blood investigations. From ECG finding and abnormal cardiac function test. Diagnosis of myocardial infarction was confirmed.

Following treatment was given

- loading dose of anti-platelet drug (Aspirin)
- loading dose of hypocholesterolemic (Statin group) drug
- Fibrinolytic drug (streptokinase)
- 50% dextrose saline with Thiamine (Vitamin B1)

After complete management and recovery after 7 days of admission in hospital, at time discharge from hospital, physician advised to take medicines regularly and to take more amount of fruit and fiber food.

- Random Blood Sugar = 30 mg %
- HbA1C = 9 %
- S. Cholesterol = 350 mg %
- S. Triglyceride = 250 mg %
- S. HDL Cholesterol = 25 mg %

Question

1. What are chronic complication of DM?
2. Why uncontrolled diabetic mellitus increase chances of atherosclerosis?
3. What is cardiac function test?
4. Which test will you prefer to do for diagnosis of myocardial infarction, if patient come after 4 day of onset of chest pain?
5. How statin reduce cholesterol level?
6. What is biochemical explanation of hypoglycemia?
7. Why physician asked to give injectable 50% Dextrose saline with Thiamine (Vitamin B1)?
8. What is role of fruits and fiber in chronic diabetes mellitus and atherosclerosis?
9. Why blood sample for blood sugar estimation is collected in fluoride containing vial?
10. What is re-perfusion injury ? And what is role of allopurinol to prevent it?
11. How will you calculate patient's LDL cholecterol?
12. What is role of fibrinolytic drugs (streptokinase) in myocardial infarction?

Clinical Case – 2

56 year male patient came in emergency with alter-conciuosness & haemetemesis . He was suffering from chronic cirrhotic liver disease due to chronic alcoholism. On examination , it was found that he has edema on both lower limb, fluid collection in peritoneal cavity (Ascites), yellowish discolouration of skin & sclera (icterus), with hypotension (decrease Blood Pressure).

On blood investigation following was found.

- Blood Glucose : 50 mg%
- Serum Protein : 5.5 gm %
- Serm Albumin : 2.0 gm%
- Serum Ammonia : Very High
- Serum Total Billirubin : 20 mg%
- APTT – Test : 60 second
- APTT – Control : 30 second
- APTT – INR : 2
- Haemogloin : 6 gm%

Ultra Sono-Graphy detected

- Cirrhosis of Liver
- Fatty Liver

Physician advise to give Following treatment

- Injection 10% Dextrose
- Injection Thiamine (B1)
- Injection Vitamin K
- Injection 10% Albumin
- Oral Neomycin (Anti-microbial, Antibiotic)
- Liq Lactulose (Laxative)
- Oral Phenylbutarate

1. Biochemical explanation about following symptoms in chronic alcoholic
 - a. Alter conciousness
 - b. Haemetemesis
2. Biochemical explanation about following signs in chronic alcoholic
 - a. Edeme
 - b. Ascites
 - c. Hypotension
3. What is hepato-renal syndrome?
4. Biochemical reason for giving following in patient of chronic alcoholic
 - a. Dextrose plus thiamine
 - b. Vitamin K
 - c. 10% Albumin
 - d. Oral Neomycin (Anti-microbial, Antibiotic)
 - e. Liq Lactulose (Laxative)
 - f. Oral Phenylbutarate

Clinical Case – 3

A 54 year old obese person come in emergency with altered consciousness level and increase respiratory rate (tachypnea) for last 4 hours.

He is having history of uncontrolled diabetes mellitus since 15 years, as he was not following any medical advice from physician. He was on insulin therapy for 3 years, but he was not taking regular dose of insulin. Patient's relative is telling that he is also having complain of weakness and decrease urine output for last 2 days.

On General examination, physician noted

- Dryness of mouth
- Pale & dry conjunctive
- Shrunken eye ball.
- Feeble (low volume) pulse
- Tachypnea (increase respiratory rate)
- Tachycardia (increase heart rate)
- Very low blood pressure (70/40 mm Hg).

Doctor makes admission in ICU and asked immediately for blood investigation.

Parameter	Value	Reference range
RBS	500 mg/dl	140 mg/dl
Serum Acetone	10 mg/dl	<1 mg/dl
Serum Creatinine	2.5 mg/dl	0.4 - 1.4 mg/dl
Blood Urea	150 mg/dl	15 - 45 mg/dl
Serum Na ⁺	120 mmol/l	135 - 145 mmol/l
Serum K ⁺	6.0 mmol/l	3.5 - 5.0 mmol/l
pH	7.1	7.35 - 7.45
pO ₂	95 mmHg	90 - 100 mmHg
pCO ₂	24 mmHg	32 - 40 mmHg
HCO ₃ ⁻ (Bicarbonate)	12 mmol/l	24 - 32 mmol/l

Diagnosed = “Diabetic ketoacidosis with acute renal failure”

Advised to following treatment.

- Inj normal saline fast I.V. (4-5 litre in 1st 24 hrs) Until systolic blood pressure reaches to normal
- Inj Human Insulin injection slow infusion I.V.As per blood sugar level
- Inj Bicarbonate 200 ml I.V.
- K⁺ Binding resin Sachets Orally.

Urinary catheterization done.

But urine output is nil

To follow below protocol for treatment of this patient.

- If RBS > 200 mg/dl ---> Give Normal Saline + Human Insulin
- If RBS < 200 mg/dl ---> Give Dextrose Saline + Human Insulin

Doctor asked to repeat following investigation during management

- RBS every 2 hourly.
- Serum K⁺ level after 4 hours.
- Arterial Blood Gas analysis after 6 hours (if require)

24 hours after admission and intensive care

- He get consciousness, normal respiration ,
- normal blood pressure & 1200 ml of urine output.
- RBS = 150 mg% with Human insulin infusion
- Serum acetone = 2 mg/dl
- Electrolyte and ABG = Normal.

He shifted to ward & remained admitted for 5 days in hospital.

On discharge, physician advises to take prescribe insulin dose regularly as well as regular follow up with FBS & PP2BS.

1. Give explanation for altered consciousness and increase respiratory rate in this case.
2. What metabolic and functional abnormality can occur due to increase acetone level?
3. Why after 24 hours serum acetone came down nearer to normal level?
4. What is patho-physiology behind decrease urine output in this patient?
5. Give comment on patient ABG report.
6. Give biochemical reason for increase K⁺ level in this case.
7. What is biochemical reason for giving dextrose saline plus human insulin infusion if RBS is below 200 mg%?
8. How bicarbonate, insulin and K⁺ binding resin reduce serum potassium level?

Medical Biochemistry – Syllabus

Medical Biochemistry encompasses any topic of biochemistry relevant to human health and diseases. As medicine is an ever expanding body of knowledge, Medical Biochemistry syllabus is continuously expanding.

At bare minimum, you are expected to get integrated knowledge of theoretical and practical aspects of following in context of the field of Medicine. In addition, newer advances in the field of medical biochemistry needs to be studied.

Carbohydrates:

Chemistry, Nutrition, Digestion, Absorption, Transport, metabolism and biochemical basis of related diseases, their treatment and prevention.
Alcohol metabolism

Amino acids and Proteins:

Chemistry, Nutrition, Digestion, Absorption, Transport, metabolism and biochemical basis of various diseases, their treatment and prevention.
Enzymes
Hemoglobin and Heme metabolism
Plasma proteins
Collagen, elastin and extracellular matrix proteins

Lipids:

Chemistry, Nutrition, Digestion, Absorption, Transport, metabolism and biochemical basis of various diseases, their treatment and prevention.
Prostaglandins

Nucleic Acids:

Chemistry, Nutrition, Digestion, Absorption, Transport, metabolism and biochemical basis of various diseases, their treatment and prevention.

Genetics

DNA and RNA structure and functions
Genome and Chromatin
Replication, Transcription, Genetic code and Translation
DNA Damage and repair
Mutations
Recombinant DNA Technology
Cell cycle and its regulation
Biochemistry of cancer
Biochemical basis of genetic diseases, their treatment and prevention.

Integration of metabolism:

Bioenergetics

Cellular Respiration

Interrelationship among metabolic pathways.

Biochemical basis of related diseases, their treatment and prevention.

Vitamins:

Chemistry, Nutrition, Digestion, Absorption, Transport, metabolism and biochemical basis of various diseases, their treatment and prevention.

Minerals:

Chemistry, Nutrition, Digestion, Absorption, Transport, metabolism and biochemical basis of various diseases, their treatment and prevention.

Water and pH:

Water biochemistry and biochemical basis of related disorders.

Blood buffers, regulation of blood pH and biochemical basis of related disorders.

Xenobiotics:

Chemistry, Metabolism and excretion of xenobiotics.

biochemical basis of related disorders

Tools for study of Biochemistry:

Colorimetry

Chromatography

Electrophoresis

ELISA

RIA

PCR and blotting techniques

Biochemistry of supramolecular structures (overlapping above topics):

Biochemical characteristics of various organelles, cells, tissues and organs e.g. Mitochondria, peroxisomes, general cell structure, RBC, Liver, Brain, Heart, Skeletal muscles etc.

Subject distribution and paper style

Paper distribution:

Paper 1:

Chemistry, digestion, absorption and metabolism of *Carbohydrate, Lipid, Water, pH, Minerals*

Paper 2:

Chemistry, digestion, absorption and metabolism of Protein(including hemoglobin, plasma proteins and enzymes), Nucleic acids including genetics, Vitamins, Xenobiotics

Note: *Overlapping common topics are acceptable in any paper e.g integration of metabolism, nutrition, tissue and organ biochemistry, biochemistry laboratory techniques, biochemistry of microorganisms (e.g HIV), environmental biochemistry and Cancer.*

Paper style(paper 1 and 2)

Section 1

Q-1 Short notes (2 out of 3)	08 marks
Q-2 Describe in brief (4 out of 6)	12 marks
Q-3 Write answer in few line(5 out of 6)	05 mqrks

Section 2

Q-4 Case with 5 questions	10 marks
Q-5 Answer in few lines(5 out of 7)	10 marks
Q-6 Write answer in few line(5 out of 6)	05 marks

Important Short Notes

General

1. Blood buffer mechanism
2. Renal Buffer mechanism for acid base balance.
3. Arterial Blood Gas Analysis & interpretation.
4. H₂O₂ – Myeloperoxidase (MPO) – Halide system of ROS (reactive Oxygen species)
5. Fluidic Model of Cell membrane
6. Type and Example of Transport mechanism.
7. Primary & Secondary cause of Hyperuricemia (Gout)
8. Chemi-osmotic hypothesis.
9. Energy production in TCA cycle.
10. Uncouplers of Oxidative phosphorylation
11. Principle, Type and utility of Electrophoresis.
12. Principle, Type and utility of ELISA.
13. Principle and utility of Colorimeter
14. Biochemical changes in Liver, Adipose tissue and muscle in well fed state
15. Biochemical changes in Liver, Adipose tissue and muscle in well fasting.

Carbohydrate

16. Mucopolysaccharide
17. Proteoglycans
18. Digestion & absorption of carbohydrate
19. Lactose intolerance
20. Energy production of Glycolysis
21. Regulation of Glycolysis
22. Amphibolic role of TCA cycle
23. Significant of NADPH
24. Significant of HMP Shunt pathway
25. Substrate & Regulation of Gluconeogenesis
26. Sorbitol Pathway
27. Polyol pathway and its significant
28. Effect of alcoholism on gluconeogenesis, oxidation of fatty acid & TCA cycle.
29. Diagnostic criteria for Diabetes Mellitus
30. Define and significant of Glycate haemoglobin
31. Metabolic alteration in Diabetes Mellitus
32. Acute and Chronic complication of Diabetes Mellitus
33. Biochemical explanation of Diabetic Ketoacidosis
34. Define C-peptide & its significant.
35. Define and significant of Glycated (HbA1c) haemoglobin
36. Advance glycated end product.
37. Type of Diabetes Mellitus . Explain LADA, MODY, Gestational diabetes & Secondary diabetes
38. Von Gierke's Disease
39. G6PD deficiency
40. Fate of Acetyl CoA.
41. Ketone body synthesis & utilization
42. Alcohol metabolism
43. Epimer & Isomer

Lipid

44. Name and Significant of Essential Fatty acid
45. Type, Difference & clinical significant of saturated & unsaturated fatty acid.
46. Fate of cholesterol
47. Ketogenesis and ketolysis
48. Metabolism of LDL
49. Regulation of LDL receptor
50. Formation of Eicosanoid & explain its inhibitor with significance.
51. Pathogenesis of atherosclerosis in context of 'oxidised LDL'.
52. Type and differentiation of oxidation of fatty acid
53. Rancidity of Fatty acid
54. Liposome & Micelle
55. Function of Phospholipids
56. Lung surfactant
57. Role of Phospholipase A2 of Snake venom in RBC lysis.
58. Role of phospholipid in signal transmission
59. Lipid digestion –absorption.
60. Significance and Regulation of Cholesterol.
61. Risk factor for Atherosclerosis
62. Prevention of Atherosclerosis
63. Type and Function Lipoproteins
64. Type and function of Apo-lipoproteins
65. Metabolism of HDL
66. Reverse cholesterol transport.
67. Role of Lipoprotein-a in Atherosclerosis
68. Energy production of long (16 carbon) saturated fatty acid through Beta oxidation
69. Carnitine shuttle
70. Assessment, Metabolic changes and influencing factors of obesity.
71. Cause of Fatty liver
72. Name the Lipotropic Factor. Explain its effect

Protein and Amino acid

73. Zwitter ion
74. Functional classification of protein
75. Protein structural –functional relationship
76. Primary, Secondary, Tertiary and Quaternary structure of Protein.
77. Define Chaperon & Prion protein
78. Protein folding & unfolding.
79. Define Protein Denaturation. Give its significant & causative factor
80. Digestion & Absorption of Protein
81. Urea cycle & Ammonia detoxification.
82. Define & give significant of transamination, transdeamination & deamination.
83. Overview of tyrosine & phenylalanine
84. Fates of Tyrosine & Phenylalanine & its related disorder
85. Biochemical explanation of Phenylketonuria
86. Biochemical explanation of Albinism & Alkaptonuria
87. Fates of Tryptophan & its related disorder.
88. Maple Syrup Urine Disease

89. Folate trap
90. Collagen-Homocystineuria-Ectopia lentis
91. Nitrogen disposal-GDH and Alpha ketoglutarate
92. Role of Glutathione & NADPH for maintain RBC membrane
93. Fates of Glycin
94. Fates of Glutamic acid
95. Transport and Detoxification of Ammonia
96. Role of 2-3 BPG on oxygen diffusion-dissociation and effect during hypoxia
97. Mechanism of the Halden & Bohr effect
98. Developmental changes in Hemoglobin gene expression from intrauterine life to adult.
99. Regulation of Hemoglobin synthesis.
100. Haemoglobin degradation pathway & it's related disorder.
101. Types , Causes and differentiation of Jaundice by serum and urine examination.
102. Haemoglobin synthesis pathway & it's related disorder.
103. Define Porphyria. Explain Causes, Clinical Feature and diagnosis of Acute intermittent porphyria and Congenital erythropoietic porphyria.
104. Molecular and Biochemical explanation for pathogenesis of Sickle cell disease
105. Molecular and Biochemical bases of Thalassemia.
106. Define and explain cause & effect of Met-haemoglobinemia
107. Transport Plasma proteins
108. Storage proteins

Enzyme

109. Write and Explain Factor affecting enzyme activity with example.
110. Explain First order & zero order enzyme kinetics.
111. Explain Difference in Function of Glucokinase and Hexokinase on bases of it's V_{max} and K_m .
112. Difference between Competitive inhibition and Noncompetitive inhibition.
113. Diagnostic importance of isoenzyme
114. Type of Enzyme Inhibition. Explain with example.
115. Regulation of Enzyme activity
116. Define Co-Enzyme, Cofactor , Apo-Enzyme , Prosthetic group & Holoenzyme
117. Enumerate Liver Function Test & Write it's significant.
118. Enumerate Cardiac Function Test & Write it's significant.

Nutrition & Vitamin

119. Protein Energy Malnutrition (PEM)
120. Difference between Kwashiorkor & Marasmus
121. Factor affecting Basal metabolic rate
122. Assessment of obesity.
123. People having Mediterranean diet show low incidence of CHD
124. Clinical significance of Dietary fibre
125. Vitamin B12 & folic acid deficiency can cause hyperhomocysteinemia
126. Folate trap
127. Effect of Warfarin & Dicoumarol on Vitamin K metabolism
128. Function of Vitamin K
129. Function of vitamin C

- 130.**Role of Active form of Vitamin B1 (Thiamine) in metabolism and give it's significance.
- 131.**Visual cycle of Vitamin A
- 132.** Name and write clinical manifestation occur in Vitamin A deficiency.
- 133.**Metabolism, Function and Clinical significance of Vitamin D
- 134.**Homoestasis changes in calcium , vitamin D & parathyroid hormone in case of renal failure.
- 135.**Regulation of calcium.
- 136.**Hypocalcaemia
- 137.**Mucosal block theory of iron absorption.
- 138.**Iron deficiency Anemia
- 139.**Type and clinical features of Beriberi.
- 140.**Pernicious anaemia.
- 141.**Metabolic changes during starvation
- 142.**Metabolic role of Vitamin B12.
- 143.**Name Riboflavin (FAD) & Niacin (NAD⁺ & NADP⁺) dependant enzymatic reaction.

Molecular

- 144.**Type and Watson & Crick Model Of DNA
- 145.**Organisation of eukaryotic DNA.
- 146.**t-RNA.
- 147.**Degeneracy & wobbling phenomena
- 148.**Genetic codon
- 149.**Molecular basis of Sickle cell anaemia.
- 150.**Type of DNA polymerase & specify it's fuction.
- 151.**Name & role of the component of the DNA replication fork
- 152.**DNA repair mechanism.
- 153.**Define Telomer & Telomerase. It's significant
- 154.**Effect and Type of Mutation with examples.
- 155.**Initiation of Transcription
- 156.**Post-transcription modification.
- 157.**Post translation modification.
- 158.**Protein synthesis inhibition by drugs.
- 159.**Salvage pathway of Purine synthesis and related disease.
- 160.**Lysch Nyhan Syndrome
- 161.**Adenosine deaminase deficiency.
- 162.**Lac operon
- 163.**Procedure & Significant of PCR
- 164.**Significant of RFLP in diagnosis of Sickle cell disease
- 165.**Microarray
- 166.**Recombinant DNA Technology
- 167.**DNA Library
- 168.**Uric acid synthesis and its inhibitors.
- 169.**Causes and management (biochemical aspect) of Gout
- 170.**Type and function of Topoisomerase
- 171.**DNA dependent RNA polymerase.
- 172.**Ribozymes
- 173.**Define Autosomal dominant & Autosomal recessive & Draw pedgree chart.